

## Product datasheet for **SC316803**

### PHACTR2 (NM\_001100164) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PHACTR2 (NM_001100164) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHACTR2
Synonyms:	C6orf56
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC316803 representing NM\_001100164.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGGCCAGACCTCGGTGTCCACGCTGTCCCGCAGCCCGGCAGCGTTGACGGACTGGACAAAGCTTCT
ATAGCAAACCTCAGATGGCCCCACAGCAGTTCACAAACACCTCCCTTCAAAGAAAAGGGAAACTATCC
ACCATTGGTAAAACTTTAAGCCTTGGAAATGGAGGAAAAAGAAGACCAGCGACAAATTTAGAGAAACC
TCGGCAGTATTAGAAAGGAAGATATCCACACGACAAAGTAGAGAGGAGCTGATAAGAAGGGGAGTGCTT
AAGGAATTGCCTGATCAAGATGGAGATGTAAAGTTAACTTTGAAAATCAAACGGGCACATGATACCC
ATCGGAGAGGAATCTACCCGAGAGGAAAATGTAGTAAAGTCTGAAGAAGGTAATGGCTCTGTATCTGAA
AAAACACCACCTCTGGAGAACAGGCAGAAGATAAGAAAGAGAACACTGAAAACCACTCTGAAACACCG
GCAGCTCTGCTCTACCTCCTTCTGCTCCTCTAAGCCTAGACCCAAACCTAAACCCAAAAAATCACCT
GTGCCTCCGAAAGGGCCACTGTGGGGCCAGCCACAAAGGTGATGAAGTGCCTCCCATTAATAAAAAAT
ACCAAGGCTCCTGGTAAGCAGGCCCCCGTCCCTCCACCAAGCCAGCAAGCCGAAACACGACCCGAGAG
GCTGCTGGCTCCTCTATTCAAAAAAACAACCTGGCTCTAAAGCATCAGCTTCGCCATCCACTTCATCC
ACCTCATCTCGTCCCAAAGCTTCAAAGGAGACAGTTTCTAGCAAAGCAGGGACAGTGGGGACCACCAAG
GGCAAGAGAAAAACTGACAAGCAGCCAATAACTTCTCACCTGTCTCAGACACAACAACCTCTGGCACA
TCCGACCTGAAAGGAGAGCCTGCAGAGACCAGAGTGGAGAGTTTCAAACCTCGAACAGACTGTCCCTGGA
GCTGAGGAGCAGAACACAGGCAAATCAAGTCCATGGTCCCTCCACCCCTGTGGCTCCAGCACCTTCT
CCTCTGGCCCCCTCTCCCTCTTGAGGATCAGTGCATTACTGCCTCAGACACTCCAGTTGTCCTCGTC
AGCGTTGGAGCTGACCTGCCCGTCTCTGCCTTAGACCCAAGTCAGCTTCTTTGGGCTGAAGAGCCGACG
AACAGAACCACCTCTACTCAGGCAGTGGCTTAAGTGTTAACAGAGAAAAAGCAAAATGTTTCACAACC
AAAGAGGAGCTGGGGAAGACAGTGCCTCAGCTACTGACTCCTGGGCTGATGGGCGAATCTTCAGAATCC
TTTAGTGCTCAGAAGATGAAGGCCACAGGAATACCAAGCCAATGACTCTGACTCGGACGGCCTATC
TTGTACACCGATGATGAGGACGAAGACGAAGATGAGGATGGCAGTGGAGAAAAGTCTTTGGCAAGTAAA
ATACGCCGGAGGGATACTCTTGCTATCAAACCTGGCAACAGACCATTAAGAAAGAACTAGAGGACAAA
AACATCTTGACAGCTACATCTGAAGAAGAGAGGCAGGAAATCCGACAACAAATTTGAAACCAAGTTAGTC
AGGAGGCTGAGCCAGAGGCCACAACCTGAAGAATTAGAGCAAAGAAACATCTAAACAAAAGAATGAA
GAAGAGGAACAAGAAGCAAAAAATGGAACCTAACGCAGACTCAGCAGAAAGCTCAGCCTGAGACCCACA
GTGGCAGAGCTACAAGCTCGAAGGATCCTGCGATTTAACGAGTATGTAGAAGTCACGGATTCCTCTGAC
TATGACCCCGCAGCAGACAAGCCCTGGGCCAGGTTAACACCTGCAGACAAGGCAGCAATAAGGAAAGAA
CTAAATGAATTTAAAAGCACAGAAATGGAAGTTCATGAAGAGAGTCGACAGTTTACAAGGTTTCATCGT
CCATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:**

**ACCN:** NM\_001100164

**Insert Size:** 1938 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001100164.1</u>
<b>RefSeq Size:</b>	9642 bp
<b>RefSeq ORF:</b>	1938 bp
<b>Locus ID:</b>	9749
<b>UniProt ID:</b>	<u>O75167</u>
<b>Cytogenetics:</b>	6q24.2
<b>MW:</b>	70.7 kDa