

Product datasheet for **SC316822**

FBXL19 (NM_001099784) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FBXL19 (NM_001099784) Human Untagged Clone
Tag:	Tag Free
Symbol:	FBXL19
Synonyms:	CXXC11; Fbl19; JHDM1C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316822 representing NM_001099784.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGGTATGAAAGTCCCCGGAAGGGGAATCTGGGCCCTCGGCGTTGCTGACGCCCCCAATGTCGTCG
AGCAGCCGGGGCCGGGGCCGGAGCGCGCCGACGCCGAACCCGCTGCCGCCGCTGCCGGCCCTGTGTG
CGAACTGAGTGCGGGGATTGCCACTTCTGCCGAGACATGAAGAAGTTCGGGGGGCCCGGGCGCATGAAG
CAGTCGTGCCTGCTCCGGCAGTGCACTGCCCCGTGCTCCACACACAGCTGTGTGCCTCTGTGTGGG
GAGGCTGGGAAGGAGGACACGGTGGAGGGAGAGGAAGAGAAATTTGGTTTGAGCCTCATGGAGTGACA
ATCTGCAACGAGATCGTCCACCCCGGCTGCCTGAAGATGGGAAGGCTGAGGGTGTCAATGCAGAG
ATCCCCAACTGCTGGGAGTGCCCTCGCTGCACCCAGGAAGGCCGACAGCAAGGATTAGGTGAGGGG
CCTGGCCGCCGTAGGGCCGACACGGCGAGGAGGGCGCCAGCTTGGGAGCGGATGGAAGCTGACAGAG
GAGCCACCGCTTCCACCGCCCCCGCCAGGCGCAAGGGCCCCCTGCTGCCGGGGCCCCCGGAGGAC
GTGCTTGGGCCCCCAACGAAAGGAAAGGGAGGCAGGGAATGAGCCTCCACCCCAAGGAAAAAGGTG
AAAGGAGGCCGAGAGAGGCACCTGAAGAAGGTGGGTGGAGACGCCTGCCTCCTCCGAGGATCGGACCCA
GGCGGCCCGGGCTGCTGCCCCCAGGGTTCTGAATCCGAGCCAGGCTTCTCATCCTGCCACCCCTGGG
CTCCCTCCCGAGAACTGGGAGAAACAAAGCCGCTTTGGCCTCTGCAGAGGGCCAGCGGTGCCGTCC
CCGTCCCCGAGAGGGAGAAGCTAGAGCGTTTCAAGCGGATGTGCCAGCTGCTGGAACGGGTGCCTGAC
ACCTCCTCTTCTCCTCGGACTCAGACTCCGACTCTTCGGGCACATCGCTGAGTGAGGACGAA
GCCCCGGCGAGGCCCGGAATGGGCGACGGCCAGCCCGGGGCGAGCTCTGGCGAGAAGGAGAACCCTGGG
GGGCGGCGGGCTGTGCGCCCTGGCAGTGGGGGGCCCTACTCAGCTGGCCCTGGGCCAGCCCCACCA
CCCCGGCTCCACAGCTGGAGCGGCACGTGGTGGCGCCCCCGCTCGAAGCCCTGAGCCCGACACTC
CCCTTGGCTGCTGGATCCGACCACCCCTGCCCGGGCCGCTGGCTTCGCGTCTTCCAGCACCTCGGG
CCGCGGGAGCTGTGTATCTGCATGCGAGTCTGCCGAATTTGAGCCGCTGGTGTATGACAAGCGTCTG
TGGCCTCGAATGGACCTGAGCCGGCGGAAGTCACTGACCCCGCCATGCTCAGTGGTGTGGTTCGCCGC
CAGCCCCGTGCCCTGGACCTCAGCTGGACAGGTGTCTCAAGAAGCAGCTCATGTGGCTTCTGAACCGA
CTACAAGCCTGCAGGAGCTGGTGTCTCTGGCTGCTCTGGCTCTCTGTCTTGGCCTGGGCTCAGCC
CCACTGCCAGCCTGCGGCTCCTGGACCTCCGCTGGATCGAGGATGTTAAAGACTCCAGCTCCGGGAG
TTGCTGCTGCCTCACCAGACACCAACCAGGGCAAACAGAGAGCCGTGGTGGCTGCAGGGGTGGCA
GAACTGCGTCTGGCAGGTTTGGAGCTGACAGATGCCTCCCTGCGTCTCCTGCTGCGTCACGACCCAG
CTGAGCGCCCTGGACCTGAGCCACTGCGCCACGTGCGGGACCCAGTGTTCACCTCCTACGGCCCCC
ACGTCCCCACTCCGCGAGACCCTGGTGCACCTCAATCTTGCTGGTTGCCACCGCCTAACGGACCACTGC
CTCCCGCTGTTCCGCGCTGCCCTCGTCTACGCCGCTAGACCTGCGCTCCTGCCCCAGCTCTACCC
GAAGCTTGTGCCCGGTGGCAGCTGCCGGGCCCTGGCCCTTCCGCTGCCCTGAGGAGAAGCTGCTT
CTCAAGGACAGCTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI
ACCN: NM_001099784
Insert Size: 2085 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001099784.2</u>
RefSeq Size:	3814 bp
RefSeq ORF:	2085 bp
Locus ID:	54620
UniProt ID:	<u>Q6PCT2</u>
Cytogenetics:	16p11.2
Protein Families:	Druggable Genome, Transcription Factors
MW:	75.7 kDa
Gene Summary:	This gene encodes a member of the Skp1-Cullin-F-box family of E3 ubiquitin ligases. The encoded protein is reported to bind to the transmembrane receptor interleukin 1 receptor-like 1 and regulate its ubiquitination and degradation. This protein has been linked to the regulation of pulmonary inflammation and psoriasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013] Transcript Variant: This variant (1) represents the longer transcript and encodes the longer isoform (1).