

## Product datasheet for **SC316823**

### LRRN3 (NM\_001099660) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LRRN3 (NM_001099660) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRN3
Synonyms:	FIGLER5; NLRR-3; NLRR3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC316823 representing NM\_001099660.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAGGACATGCCACTCCGAATTCATGTGCTACTTGGCCTAGCTATCACTACACTAGTACAAGCTGTA
GATAAAAAAGTGGATTGTCCACGGTTATGTACGTGTGAAATCAGGCCTTGGTTTACACCCAGATCCATT
TATATGGAAGCATCTACAGTGGATTGTAATGATTTAGGTCTTTAACTTTCCAGCCAGATTGCCAGCT
AACACACAGATTCTTCTCTACAGACTAACAATATTGCAAAAATTGAATACTCCACAGACTTTCCAGTA
AACCTTACTGGCCTGGATTTATCTCAAAACAATTTATCTTCAGTCACCAATATTAATGTA AAAAAGATG
CCTCAGCTCCTTTCTGTGTACTAGAGGAAAAACAACTTACTGAACTGCCTGAAAAATGTCTGTCCGAA
CTGAGCAACTTACAAGAACTCTATTAATCACAACCTGCTTTCTACAATTTACCTGGAGCCTTTATT
GGCCTACATAATCTTCTCGACTTCATCTCAATCAAATAGATTGCAGATGATCAACAGTAAAGTGGTTT
GATGCTCTTCAAATCTAGAGATTCTGATGATTGGGAAAAATCCAATTATCAGAATCAAAGACATGAAC
TTTAAGCCTCTTCAATCTTCCGACGCTGGTTATAGCTGGTATAAACCTCACAGAAATACCAGATAAC
GCCTTGGTTGGACTGGAAAACTTAGAAAGCATCTCTTTTACGATAACAGGCTTATTAAAGTACCCCAT
GTTGCTCTTCAAAAAGTTGTAATCTCAAATTTTGGATCTAAATAAAAAATCCTATTAATAGAATACGA
AGGGGTGATTTTAGCAATATGCTACACTTAAAAGAGTTGGGGATAAATAATGCCTGAGCTGATTTCC
ATCGATAGTCTTGCTGTGGATAACCTGCCAGATTTAAGAAAAATAGAAGCTACTAACAACCTTAGATTG
TCTTACATTCACCCCAATGCATTTTTAGACTCCCAAGCTGGAATCACTCATGTGAACAGCAATGCT
CTCAGTGCCCTGTACCATGGTACCATTGAGTCTCTGCCAACTCAAGGAAATCAGCATAACAGTAAC
CCCACAGGTGTGACTGTGTATCCGTTGGATGAACATGAACAAAACCAATTCGATTCATGGAGCCA
GATTCAGTGTTTTGCCTGGACCCACCTGAATCCAAGGTGAGAATGTTTCGGCAAGTGCATTTAGGGGAC
ATGATGGAAATTTGTCTCCCTCTTATAGCTCCTGAGAGCTTTCCTTCTAAATCTAAATGTAGAAGCTGGG
AGCTATGTTTCCTTTCACTGTAGAGCTACTGCAGAACCACAGCCTGAAATCTACTGGATAACACCTTCT
GGTCAAAAACCTTGCCTAATACCCTGACAGACAAGTTCTATGTCCATTCTGAGGGAACACTAGATATA
AATGGCGTAACTCCCAAAGAAGGGGTTTATATACTTGTATAGCAACTAACCTAGTTGGCGCTGACTTG
AAGTCTGTTATGATCAAAGTGGATGGATCTTTCCACAAGATAACAATGGCTCTTTGAATATTAATA
AGAGATATTCAGGCAATTCAGTTTTGGTGTCTGAAAGCAAGTTCTAAAATCTCAAATCTAGTGT
AAATGGACAGCCTTTGTCAAGACTGAAAATTCATGCTGCGCAAAGTGTCTGAATACCATCTGATGTC
AAGGTATATAATCTTACTCATCTGAATCCATCACTGAGTATAAAATTTGATTGATATCCACCATC
TATCAGAAAAACAGAAAAAATGTGAAATGTCACCACCAAGGTTTGCACCCTGATCAAAAAGAGTAT
GAAAAGAATAATACCACAACACTTATGGCCTGTCTTGGAGGCCTTCTGGGGATTATTGGTGTGATATGT
CTTATCAGCTGCCTCTCTCCAGAAATGAACTGTGATGGTGGACACAGCTATGTGAGGAATTACTTACAG
AAACCAACCTTTGCATTAGGTGAGCTTTATCCTCCTCTGATAAATCTCTGGGAAGCAGGAAAAAGAAAA
AGTACATCACTGAAAGTAAAAGCAACTGTTATAGGTTTACCAACAATATGTCCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

**Restriction Sites:** SgfI-MluI

**Plasmid Map:** □

**ACCN:** NM\_001099660

**Insert Size:** 2127 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001099660.1</a></u>
<b>RefSeq Size:</b>	3744 bp
<b>RefSeq ORF:</b>	2127 bp
<b>Locus ID:</b>	54674
<b>UniProt ID:</b>	<u><a href="#">Q9H3W5</a></u>
<b>Cytogenetics:</b>	7q31.1
<b>Protein Families:</b>	Transmembrane
<b>MW:</b>	79.4 kDa