

Product datasheet for **SC316827**

MROH6 (NM_001100878) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MROH6 (NM_001100878) Human Untagged Clone
Tag:	Tag Free
Symbol:	MROH6
Synonyms:	C8orf73
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316827 representing NM_001100878.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGGCTGGGGGTGTGTGGGGCCGAGCCGGGCCGGGAGGCTCCCGTGGGGGCTCTAACCTGACAGCA
CTGACTGAAGGAATCCGGGCCAGGCAGGGGCAGCCCCAGGACCCCTTCCGAGGCCCTCAGCCCAAG
TCCTGGGAGGTCAAACCTGAGGCTGAGCCACAGACCCAGGCACTACCGCCCCCTCTGAGGCAGAGCCT
GGACGTGGGGCCACCGTCCCTGAAGCTGGCAGCGAGCCCTGCTCCCTCAACAGTGCCTGGAACAGCC
CCTGAGGGGCCCCACCAGGTTCCCCAGAGTTCTTGGGAGGAGGAGTTCTTGCCGACCTCGCGTTGTAC
ACGGCTGCCTGCCTGGAGGAGGCTGGCTTTCAGGGACCCAGGCGACAGTGTACCCCTGTCTCAGCC
CTGGAGGGCCGGGGCAGCGGTTGGAGGACCAGGTGCATGCTCTGGTGCCTGGGCTGTGGCGCAGGTG
CCCAGCCTAGCGAGGGGAGGCCCTGGAGGGCGGCCCTGCGAGTGTGAGCGCACTGGCCCTGGAGCAT
GCGCGGGACGTGGTGTGTGCGCTGCTACCCGCTCTCTGCCCGCCGATCGGGTAGCAGCCGAGCTCTGG
CGCAGCCTAAGCCGTAACCAGCGTGTAAATGGGAGGTGTGGTGAACCTGTGTGGGCGCTGAAGGT
GCTTCGGGGCCGAGCCCCAGGCACTGGCGGCCACACGTGCTCTTGGGGAGATGTGGCTGTTTCGGGC
TGCGTGGGAGCCACGAGGGGCTTCTACCCACATCTGCTGCTTGCGCTGGTCAACAGCTGCACAAGCTG
GCCCCAGCCCCGTGCTCCCCGACATGCCAAGATTTGGGTTCTGTCCCACCGAGGGCCACCACATAGC
CATGCCAGCTGTGCTGTGGAGGCCTTGAAGGCGCTGTACCCGGGATGGAGGCCGATGGTGGTCACG
TGCATGGAGCAGGCAGGAGGCTGGAGGAGGCTGGTGGGAGCCACACCCACCTGGAGGGCGTCTGCTG
CTGGCCAGTGCTATGGTGGCACATGCCGACCACACCTGCGAGGCCTCTTCGAGACTTGCTCCCTCGG
CTTCGACGCGGACGACCCGACGCTCTACGGCTATGGCTTCTTACAGGGCTGTTGAGAGCCGG
CCCACCGCAGGCTCCTGCGGGAGGAGGTATCCTGGAGCGACTCCTACCTGGCAGGAGACCCGAA
CCCACTGTGCGCTGGTTGGGCTGCTGGGCTGGGCCACCTCGCGCTGAATCGCAGGAAGGTGCGGCAC
GTGAGCACGCTGCTGCCGCGCTCCTGGGCGCACTGGGCGAAGGCGACGCGCGGCTCGTGGGTGCAGCG
CTGGGCGCCCTGAGGAGGCTCCTGCTGCGGCCCGGGCGCTGTGCGGCTCCTGAGCGCGGAGCTGGGA
CCGCGCCTCCCTCCGCTACTGGACGACACACGGGACTCAATCCGCGCTCGGCCGTGGGCTCCTGGG
ACTCTGGTGGCGGGGGCGGGGCGGCTCCGGCTGGGCTCCGCGGCCCTGCGGAAGCTGGTGTG
CAGAGTCTCGTCCGCTGCTGTCGCTGCATGACCCAGCAGGAGCGTGTGAGAGCTCAGAGTGG
ACCCTGGCCCGTGTGACCACGCTTTTGTGGGCTGCTGGAGGAGTTGGTACCCTGGCCCACTAT
GACAGCCCCGAGGCCCTGAGCCACCTGCTGCCGCTGGTTACCGGATACCCAGGCCACGTGCCAAC
TTCTGAGCCAGACCCAGGGCTACCTGCGGAGTCCACAGGACCCCTGCGCCGGGAGCCGCCGTGCTT
ATAGGCTTCTTGTCCACCACGCCAGCCCCGGCTGTGTCAACCAGGACCTGCTGGACTCCCTGTTCCAG
GACCTAGGGCGACTGCAGAGCGACCCCAAGCCGGCTGTGGCCGCGGCAGCGACGTGTCCGCTCAGCAG
GTGGCGATGTGGCCCGTGCCGGGGCTGCCCGCGGGCCCCGCTTCTCCGCATCGCCCCGCGCCCC
GCCCGGCCCCACCAGTCTTCCCGACAGCCCCCTTCCAGCGCGGAGCGTTCGCGGCGCTGGGGCTGC
TCCGGACCCCGCGAGCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001100878

Insert Size: 2160 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100878.1</u>
RefSeq Size:	3251 bp
RefSeq ORF:	2160 bp
Locus ID:	642475
UniProt ID:	<u>A6NGR9</u>
Cytogenetics:	8q24.3
MW:	77.2 kDa