

Product datasheet for SC316925

AADACL3 (NM_001103170) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	AADACL3 (NM_001103170) Human Untagged Clone
Tag:	Tag Free
Symbol:	AADACL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC316925 representing NM_001103170. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGGGACCTGGCCCTGATCTTCTCGCAGCAGCCTGCGTGTCTCACTAGGGGTCACTCTGTGGGTC
ATTTGCAGCCATTTTTCTACTGTGCACATCCCTGCAGCGGTTGGCCACCCTGTGAACTGAGAGTCTC
CATTGCATCTTCCAGTCTGTTGACTTGGGGATGATATTTGAGAAGCTCAGAATCTGTTCTATGCC
CAATTTTTCTGTTTCATGCAAGATCTGCCTCCGCTAAAGTATGACCCCGATGTTGTGGTCACGGATTT
CGCTTTGGGACAATCCCTGTGAAGCTGTACCAACCAAGGCATCCACCTGACCCCTGAAGCCTGGCATC
GTGACTACCACGGTGGCGGGGGCGTCATGGGGAGTTTGAAACCCACCATGGCATATGCTCTCGTTTG
TGCAAGGAGAGTGACTCCGTGGTTCTGGCAGTTGGTTACCGCAAGTTACCTAAGCATAAGTTTCCAGTG
CCAGTAAGAGACTGCTTGGTGGCCACCATCCACTTCTGAAGTCCTGGATGCATATGGAGTGGATCCA
GCCCGGTTGTGGTCTGCGGTGACAGTTTCGGAGGGGCAATAGCCGAGTGGTTTGTCAACAACCTGTG
GACAGGCCAGATCTGCCCGGATCCGGGCTCAGATCCTGATCTATGCCATTCTCCAAGCCCTGGATTTA
CAAACCCCTTCGTTTCAACAGAGGAAAAACATCCCACTGCTCACCTGGAGTTTCTCTGCTACTTTTTT
TTTCAAACCTGGATTTAGCTCCTCCTGGCAAGAGGTCATCATGAAAGGTGCCATTTGCCTGCTGAA
GTCTGGGAAAAGTACAGAAAGTGGTTGGGCCAGAAAACATCCCTGAGAGGTTTAAAGGAGGGGTTAC
CAACTGAAGCCCATGAGCCATGAATGAAGCTGCTTACTTGAAGTAAGTGTTCCTGGATGTGATG
TGCTCGCCCTGATTGCAGAAGTACATAGTGTCTCAGCTCCCGAAACCTGCATCGTGAGCTGTGAG
TATGATGCTCTCCGGGACAATCACTGTTGTACAAGAAAAGGCTGGAAGACCTGGGAGTGCCCGTGACC
TGGCACCATATGGAGGATGGTTTCCATGGAGTGCTCAGGACCATGACATGAGCTTCTTGCACTTTCCC
TGCTCCATGAGAATCTGAGTGCATTAGTTCAATTTGTAAGGGGACTGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: SgfI-MluI

Plasmid Map:



[View online »](#)

ACCN:	NM_001103170
Insert Size:	1224 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001103170.2</u>
RefSeq Size:	4049 bp
RefSeq ORF:	1224 bp
Locus ID:	126767
UniProt ID:	<u>Q5VUY0</u>
Cytogenetics:	1p36.21
MW:	46.2 kDa