

Product datasheet for **SC316954**

SNX18 (NM_001102575) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	SNX18 (NM_001102575) Human Untagged Clone
Tag:	Tag Free
Symbol:	SNX18
Synonyms:	SH3PX2; SH3PXD3B; SNAG1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316954 representing NM_001102575.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGATCGCC
ATGGCGCTGCGCGCCCGGGCGCTGTACGACTTCAGGTCGGAGAACCAGGAGAGATCTCGTGCAGAG
CACGAGGTGCTGAGCCTGTGCAGCGAGCAGGACATCGAGGGCTGGCTCGAGGGGGTCAACAGCCGCGC
GACCGCGCCTCTTCCCGCCTCTATGTGCAGGTGATCCGCGCCCGAGCCTGGCCCGCGGGAGAC
GGCGGCCCGGCCCGCCCGCCGCTACGCCAATGTGCCCGCGGGGCTTCGAGCCCTGCCTGTGCGC
CCCCCGCCTCTTCAAGCCGCCGCTGACGCCTCCAGGCGCTGCTGCAGCCACAGCAGGCGCCGCT
CCGAGCACCTTCCAGCCGCCGGCGGGCTTCCCGTACGGCGGGGGCGCCCTGCAGCCGTCGCTCAG
CAGCTCTACGGCGGTACCAGGCCAGCCAAGGCAGCGATGATGACTGGGACGACGAGTGGACGACAGC
TCCACGGTGGCGGACGAGCCGGCGCTCTGGGACGCGGAGCATACCCGGACCTCGACGGCTCGTCTCG
GCGGGTGTGGCGCAGCCGGCCGCTACCGCCTGTCCACGCGCTCCGACCTGTCCGTTGGTTCGCGCGC
GGCTCGGTCCCCCGCAGCACCCCGTGGGGCCCAAGAGCTCGGCCACCGTGAGCCGCAACCTCAAT
CGTTCTCCACCTTCGTCAAGTCCGGCGGGGAGGCCTTCGTGCTGGGGGAGGCGTCAGGCTTCGTGAAG
GACGGGGACAAGCTGTGGTGGTGTGGGGCCCTATGGCCCCGAGTGGCAGGAGAACCCTACCCGTTT
CAGTGCACCATCGACGACCCACCAAGCAGACCAAGTTCAAGGGCATGAAGAGCTACATCTCCTACAAG
CTGGTGCCACGCACACGAGGTGCCGGTGCATCGGCGCTACAAGCACTTCGACTGGGTGTACGCGCGC
CTGGCGGAGAAGTTCGCGGTCACTCCGTGCCCCACCTGCCCGAGAAGCAGGCCACCGGCCGCTTCGAG
GAGGACTTCATCTAAGCGCAGGAAGGGCCTGATCTGGTGGATGAACCACATGGCCAGCCACCCAGTG
CTGGCGCAGTGCAGCTTCCAGCACTTCCTGACGTGCCCGCAGCAGCACCGACGAGAAAGCCTGGAAG
CAGGGCAAGAGGAAGCCGAGAAGGACGAGATGGTGGGCGCAACTTCTTCTGACCTTAGCACGCCCC
CCGCGCGCTGCCCTTGACCTGCAGGAGGTGGAGAGCAAGATCGACGGCTTCAAGTGCTTACCAAGAAG
ATGGACGACAGCGCGCTGCAGCTCAACCACAGGCCAACGAGTTCGCGCGCAAGCAGGTGACCGGCTTC
AAAAAGGAGTATCAGAAGGTGGGCCAGTCTTCCGCGGCCTCAGCCAGGCCTTTGAGCTGGACCAGCAG
GCCTTCTCGTGGGCTGAACCAGGCTATCGCCTTACCAGGAGATGCCTATGACGCCATTGGCGAGCTC
TTCGCGGAGCAGCCAGGACGACTGGATCCCGTTCATGGACCTATTAGCGCTGTATCAGGGCATCTG
GCTAATTCCCGACATCATCCAGTTCAGAAAGGAGCTTTACCAAAGTCAAGGAGAGTAGGCGACAC
GTGGAGGAAGGAAGATGGAGGTGCAGAAGGCTGACGGCATTGAGATCGCTGTAACACTATTTCTTT
GCCACTTTGGCTGAAATTCACCACTTCCATCAAATTCGAGTGAGAGACTTAAATCACAGATGCAGCAT
TTCTTACAACAACAATAATATTTTTCCAAAAGTTACCCAGAAGTTGGAAGAAGCTCTTCAAAATAT
GATAGTGTTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map:

ACCN: NM_001102575

Insert Size: 1875 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001102575.1</u>
RefSeq Size:	5238 bp
RefSeq ORF:	1875 bp
Locus ID:	112574
UniProt ID:	<u>Q96RF0</u>
Cytogenetics:	5q11.2
MW:	69 kDa
Gene Summary:	<p>This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a SH3 domain. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]</p> <p>Transcript Variant: This variant (1) represents the longer transcript but encodes the shorter isoform (a).</p>