

## Product datasheet for **SC316964**

### Zfp219 (ZNF219) (NM\_001101672) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp219 (ZNF219) (NM_001101672) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF219
Synonyms:	ZFP219
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**Fully Sequenced ORF:** >SC316964 representing NM\_001101672.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGGAGGGCTCAGTCCCCGCGCCCCGAGCGGCCACTTAGCGCCGTCGCGCCGGCTTTCGACGGCGAG
TGGACTGCAGCGATACTCCAACGGGCCAGCCGTGAGCGCAGGGTCGCTCGGGATGGGAGCCGTGAGC
TGGTCTGAGAGTCGTGCAGGCGAACGGCGCTTCCCTGCCCTGTATGCGGGAAGCGCTTCCGCTTCAAC
TCTATCCTTGCTTTGCACCTGCGGGCGCACCCAGGAGCCAGGCCTTCCAGTGCCCTCACTGCGGCCAC
CGCGCGCGCAGCGGGCTCTGCTGCGCTCGCACCTGCGCACACACCAGCCGAGCGCCACGTAGTCCT
GCTGCACGCTGTTGCTGGAGTTGGAAGAGCGCGCTACTACGCGAGGCCCGACTGGGGAGAGCCCGA
AGCTCAGGGGGCATGCAGGCCACCCCTGCCACTGAGGGTCTGGCGCGGCCCCAGGCTCCTTCATCGTCC
GCCTTCCGTTGCCCTACTGCAAAGGCAAGTTTCGCACCTCGGCGGAGCGCAACGCCACCTGCACATC
CTGCATAGGCCCTGGAAGTGCGGCTGTGCAGTTTCGGCTCCAGCCAGGAGGAGGAGCTGCTGCACCAC
AGCCTGACGGCCACGGGGCTCCCGAGCGTCCCCTGGCGGCCACCTCCGCTGCGCCTCCGCTCAGCCT
CAGCCTCAGCCTCCACCCAGCCGAACCCAGATCAGTCCCCAGCCGGAGCCGGAGCCGGAGCCCGAA
CGTGAGGCAACCCGACCCAGCTCCTGCCGCTCCCGAGGAGCCCCAGCGCCTCCGGAGTTCGGCTGC
CAAGTGTGCGGCCAGAGCTTTACACAGTCTTGGTTTCTCAAGGGCCACATGCGTAAGCACAAAGCCCTC
TTCGATCATGCGTGTCCGGTGTGCGGCCGCTGCTTCAAGGAGCCCTGGTTCCTTAAGAACCACATGAAG
GTGCACGCCAGCAAGCTGGGCCACTGCGTGCCCCGGGGCTGCCTCCGGGCTGCCCGCGCCCCCAG
CCTCCTGACCTCGGCCTGCTGGCCTATGAGCCGTTGGGCCAGCGCTCCTTTGGCCCCGGCTCCACC
CCGGCCGAGCGCGTGCAGCCCCGAGCCTTTGGGCTACCTGAGCCTGCGAGCTGGCGAGGGCCGGCCC
AACGGCGAGGGTCTGAGCCCCGCTCCCGGCCGACGCTTCGGAGGCTTCGCCCCGCTGCTCTGCTC
CCGGCCCGGGCTCGCCGCAACCGTGCAGGAGCCTGAGGAAGAAGAGGAGGTGGTGGAGCCGAGGAG
GAAACCTGGGCCCGGGCAGGTGCTGGGCTCTTGGCTTCCCTGCACCCGCGCCGGGTGAGGGACCG
GGGCACTCTGCCTCTGCTGCTGGGGCCAGGCAAGATCGACCCGACGCAAGAGAATGGGCTGTTG
GTTGGAGGACCCGGCTGAAGGGGGCCGGGCCACCCGCAAGGATTGTCTTTCTGCGGAAAATCT
TTCCGCTCAGCACATCACCTCAAAGTGATCTGCGAGTGCACACAGGCGAGCGGCCCTACAAGTGTCCG
CACTGCGACTACGCGGCCACCCAGTCCGGCTCGCTCAAGTATCACCTACAGCGCCACCACCGGGAGCAG
AGGAGCGGGGCCGCCCGGGCCACCCCGGAGCCACCGCCTCCTTCCAGCGGGGTTCCGCCCCGCAA
TCTGGAGCAAGCCGTCTCCGCAGCCTGCGACCTGGGTGGAGGGCGCCTCAAGTCCCCGGCTCCTTCT
AGCGGTGCTGGGCGGGGTCCTGTCGAAAGCCCGCAGCCCTGGGAGGACCTGCACAACGGGCGAGGC
GGTGAGGCCGAACCCCTGGACCTGCTTTCGGGGCAGGGCCGGGAGGCGAGGCCGGGCTGGGGGTGCC
CTCCACCGCTGCCTTCTGCCGTTTCGCCACTGGAGCCCAGAGCTCATGGCCTTGACCTTCAAGTG
CACCACAGCCCGGGCTAGGGGCCCGGCCACCCAGGCTGACGCGTCCCCGCCCTATGCCCGAGTA
CCATCAGGAGAGACCCCTCCAGTCTTCGCAGGAAGGGGAGGAGGGCTCCGGCTGTCCAGACCCGGA
GAGGCAGGGCTGGGGGGCAAGAACGGTAG
AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGAT
ATCCTGGATTACAAGGATGACGACGATAAGGTTTAA
  
```

**Restriction Sites:** SgfI-RsrII

**Plasmid Map:** □

**ACCN:** NM\_001101672

**Insert Size:** 2169 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001101672.1</a></u>
<b>RefSeq Size:</b>	2964 bp
<b>RefSeq ORF:</b>	2169 bp
<b>Locus ID:</b>	51222
<b>UniProt ID:</b>	<u><a href="#">Q9P2Y4</a></u>
<b>Cytogenetics:</b>	14q11.2
<b>Protein Families:</b>	Transcription Factors
<b>MW:</b>	76.9 kDa
<b>Gene Summary:</b>	<p>This gene is a member of the Kruppel-like zinc finger gene family. The encoded protein functions as a transcriptional repressor of the high mobility group nucleosome binding domain 1 protein, which is associated with transcriptionally active chromatin. [provided by RefSeq, Apr 2017]</p> <p>Transcript Variant: This variant (2) differs in the 5' UTR compared to variant 1. Variants 1, 2 and 3 encode the same protein.</p>