

Product datasheet for **SC316972**

ZNF836 (NM_001102657) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF836 (NM_001102657) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF836
Mammalian Cell Selection:	Neomycin
Vector:	<u>PCMV6-Neo</u>
E. coli Selection:	Ampicillin (100 ug/mL)



[View online »](#)

Fully Sequenced ORF: >NCBI ORF sequence for NM_001102657, the custom clone sequence may differ by one or more nucleotides

```

ATGGCTCTTACACAGGGACCTTTGACATTCAGGGATGTAGCCATAGAATTCTCTCAGGAG
GAGTGGAAATCCCTGGACCCTGTGCAGAAAAGCTTTGACTGGGATGTGATGTTGGAGAAC
TACAGGAACCTGGTCTTCCTGGGTATCCTTCTAAATGTATGACCAAGGAATTACCACCA
ATAGGGAACAGTAATACAGGAGAAAAATGCCAAACAGTGACGCTGGAAAGACATGAATGT
TATGATGTTGAAAAATTTTTACTTAAGGGAAATCCAGAAAAATCTACAGGACCTTGAGTTT
CAATGGAAAGATGGTGAATAAATTATAAAGAAGTGCCAATGACCTATAAAAAACAATCTT
AATGGTAAAAGAGGTCAACATAGTCAAGAGGATGTAGAAAACAAATGTATTGAAAAATCAG
CTTACATTAAGCTTTTCAGTCACTGCTGACTGAACTGCAGAAAATTTCAAATGAAGGGAAA
ATTTATGAGTGTAAACCAATCTGAGAAGACAGTTAATAACAGTTCCCTAGTTTCACCACTT
CAAAGAATTTCTCTAGTGTCCAAACCAACATTTCTAAAAAATATGAGAATGAGTTTTTG
CAGCTGTCATTACCCACCAACTTGAGAAAACACACATTAGGGAAAAACCTTATATGTGT
AAAGGGTGTGGCAAAGCCTTTAGAGTGTCTTCAAGTCTTATTAACCATCAGATGGTACAT
ACTACAGAGAAACCTTACAAATGCAATGAATGTGGCAAAGCCTTTCATCGGGGCTCACTA
CTAACTATACATCAGATAGTCCATACAAGGGGAAGCCATATCAATGTGGTGTATGTGGC
AAGATCTTCAGACAAAATTCAGATCTTGTAAATCACCGGAGAAGTCACACTGGAGAGAAA
CCGTACAAATGTAATGAATGTGGCAAGTCTTTAGTCAAAGTTATAACCTTGCAATACAT
CAGAGAATTCACACTGGAGAGAAACCTTACAAATGTAATGAGTGTGGGAAAACCTTCAA
CAAGGCTCATGCCTCACTACACATCAGATAATCCATACAGGAGAGAAACCATATCAATGT
GATATATGTGGCAAGGTCTTCAGGCAGAATTTCTAATCTTGTAATCACAGAGAAATCCAC
ACTGGAGAGAAACCATACAAATGCAACATATGTGGAAAGTCTTTAGTCAAAGTTCCAAC
CTGGCAACTCATCAGACAGTTCATAGTGGAAACAAACCTTACAAATGTGATGAGTGTGGC
AAAACCTTTAAACGGAGCTCCAGCCTCACTACACATCAGATAATCCATACAGGAGAGAAA
CCATATACATGTGATGTATGCGACAAGGTCTTCAGTCAACGTTCCAACTTGCAAGGCAC
CAGAGAAGTCATACTGGAGAGAAACCTTACAAATGCAATGAATGTGGCAAGGTCTTCAGT
CAGACTTCACATCTTGTGGGCATCGGAGAATTCATACTGGAGAGAAACCTTACAAATGT
GATAAATGTGGTAAAGCCTTTAAACAGGGCTCATTACTCACTCGACATAAGATAATTCAT
ACCAGAGAGAAACGTTACCAATGCGGTGAATGTGGAAAGGTCTTTAGTAAAAATTCATGC
CTTGTAAGACATTTAAGAATTCATACTGGGAGCAACCTTACAAATGTAATGTGTGTGGC
AAGGTCTTCAATTACAGTGGAAACCTTTCAATTCATAAGAGAATACATACGGGAGAGAAA
CCTTTCCAATGTAATGAATGTGGCACAGTCTTCAGGAACACTCATGCCTAGCACGTCAT
CTAAGAATTCATACTGGGCAGAAACCTTACAAATGTAATGTGTGTGGCAAGGTCTTCAAT
GACAGTGGAAACCTTTCAATCATAAGAGAATTCATACTGGAGAGAAGCCGTTTCAATGT
AACGAATGCGGCAAGGTTTTAGTACTACTCATGCCTAGCACGTCATCGGAAAATTCAT
ACCGGAGAGAAACCTTACAAATGTAATGATTGTGGCAAAGCCTATACTCAGCGTTCAAGC
CTCACTAAACATCTGATAATTCATACTGGAGAGAAACCTTACAATGTAATGAGTTTGA
GGGGCATTATCCAAAGTTCAAACTTGCAAGATATCACAGAAATCCTACTGGGGAGAAA
CCACACAAATGTAGCCATTGTGGTAGAAGTCTTGTGATATAACAGGCCTGACGTACCAT
CAGAGAAGGCATACTGGAGAGATGCCATACAAATGTATTGAATGTGGCCAGGTCTTTAAT
TCCACTTCGAACCTTGCAAGGCATCGGAGAATTCACACTGGAGAGAAACCTTACAAATGT
AATGAATGTGGCAAGGTCTTTCGTATCAATCAACACTAGCACGTCATCGGAGTATTCAT
ACTGGAGAGAAACCTTACGTGTGTAATGAGTGTGGCAAAGCCTTTAGAGTGCCTTCAAT
CTGTTAATCATCAGAAAATGCATACTGGGGACAAACCTTACAAATGTAATGAATGTGGT
AAAGCTTTTATTGAAAGGTCAAAGCTGGTGTACCATCAAAGAAATCACACTGGAGAGAAG
CCATACAAATGTATTGAATGTGGCAAGGCCTTTGGCGGTTTTCTGCCTCAACAAACAC
CAAATGATTCATTCTGGAGAAAAACCTTATAAATGTAATGAGTGTGGCAAATCTTTCATT
AGTCGCTCAGGCCTCACTAAACATCAGACAAAACATACTGCAGAGAGTCTTAAAACAAAA
TTCAATGTGGAAAAGCCTTTAGATGTTCTCCTAACTTCTGGGTTCAA

```

Restriction Sites: Please inquire

ACCN: NM_001102657

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001102657.1](#), [NP_001096127.1](#)

RefSeq Size: 3333 bp

RefSeq ORF: 2811 bp

Locus ID: 162962

UniProt ID: [Q6ZNA1](#)

Cytogenetics: 19q13.41

Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]