

## Product datasheet for **SC316973**

### **NCOA7 (NM\_181782) Human Untagged Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	NCOA7 (NM_181782) Human Untagged Clone
Tag:	Tag Free
Symbol:	NCOA7
Synonyms:	dj187j11.3; ERAP140; ESNA1; Nbla00052; Nbla10993; NCOA7-AS; TLDC4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC316973 representing NM\_181782.  
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGATACCAAGGAAGAGAAGAAGGAACGGAAACAAAGTTATTTTGTCTGACTGAAAAAGAAAAACAA
GCCAAACAAAATGCAGAGACAGCCTCAGCTGTAGCTACAAGACTCATACTGGGAAGGAAGATAATAAT
ACAGTAGTTTTAGAGCCAGACAAGTGCAACATTGCTGTGGAAGAGGAATATATGACTGATGAGAAAAA
AAGAGAAAAAGTAAATCAGTTAAAGGAGATCAGGCGTACAGAACTAAAGAGATATTATAGTATTGATGAC
AATCAAAACAAAACACATGATAAAAAAGAGAAGAGATGGTGGTTCAGAAAGCCCATGGGACTATGGAA
TACTACTGCTGAAACCAGGACACCCTAACTCCATAGCACTGAAATTTAACATCACTCCCAATAAATTG
GTGGAAGTGAATAAACTTTTACACATACTATTGTTCCAGGCCAGGTCTTTTTGTGCCAGATGCCAAC
TCTCCTCCAGTACCTTAAGGCTATCATCATCCAGTCTGGTGTACTGTCTCTCCTTCATCATCAGAT
GCAGAATATGATAAATTGCCTGATGCTGACTTAGCACGAAAGGCCTTAAAACCCATTGAAAGAGTCTTA
TCGTCTACTTCTGAAGAAGATGAGCCAGGTGTGGTGAATTTTTAAAAATGAATTGTCGATACTTCACC
GATGGAAGGGTGTGGTGGCGGTGTTATGATAGTGACTCCTAACAAACATCATGTTTGACCCATATAA
TCTGATCCTCTGGTTATTGAAAATGGGTGTGAGGAGTATGGTCTCATCTGCCCATGGAAGAGGTTGTT
TCCATTGCGCTCTACAATGACATTTCTCACATGAAGATCAAAGATGCCTTGCCATCTGACCTACCTCAG
GATCTTTGTCTCTGTACAGGCCGGAGAATGGGAAGACCTGGCTTCAGAAAAGGATATCAACCCATTC
AGTAAGTTCAAATCTATCAACAAGGAAAAACGACAGCAGAATGGAGAGAAAAATTGACTTCGGATTCC
AGACCAATAGTACCTTTGGAGAAGTCCACAGGACATACACCTACAAGCCCTCAGGCAGCTCTGTGTCA
GAGAAATTAAGAAACTGGACTCCTCTAGGGAGACATCCCATGGTCTCCACAGTGACTAAGCTCAGC
AAGGAACCTTCCGACACTTCTTCTGCAATTTGAATCTACAGCCAAAGAAAACCTTTAGGGGAAGATGAT
GATTTTTGTTGACTTGAAGAAGCTTTCTTCTCAAAGTGGTGGTGAATGCACAAAAAGACACCTTGAAG
GAGTGCCTTTCTCTTGACCCAGAGGAACGAAAGAAAGCTGAGTCAAAATAAACAATTCTGCCGTGGAA
ATGCGAGTGCAGTCAAGCCATAGCCTTTTTGGGAACAGAGAATGATGTTGAACTGAAGGGGGCGTAGAT
TTAGAAACCTGTGAGAAGCAAGATATAATGCCAGAAGTGGACAAGCAGTCTGGTTCGCCAGAAAGCCGA
GTAGAAAAACACTGAACATACATGAAGATTTAGATAAAGTTAACTCATTGAATATTACCTGACTAAG
AACAAAGAAGGGCCACAGGTATCTGAAAATTTGCAGAAAACAGAAATTAAGTATGAAAAAGTATTGAA
CCAGGGGAATAGACATTACCCTTAGTAGTTCTCTTTCCAGGCGGGTATCCATAACTGAGGGCAAT
AAAGAGCCAGATAAGACCTGGGTGAAAAAGGGAGACCCCTCCCGTAAAACCTGAACTCTTCTACAGAA
GCAAAATGTGATTAAGAGGCTCTAGACTCCTTTGGAATCTACTCTGGACAACAGCTGTCAAGGTGCA
CAAATGGATAATAAATCTGAAGTTCAGTTGTGGCTGTTAAAGAGAATTCAGGTACCCATTGAAGATA
CTTCTTCAAAAAGAAAAAAGCAAGACCCACCCATGTTCTGTGCATCAAAGTGGGAAAACCAATG
AGAAAAATCCTTTGCCACTCACACTGCAGCCATGGTCCAGCAGTACGGCAAACGGAGAAAGCAGCCAGAG
TACTGGTTTGTCTTCTCGGGAGAGGGTGGATCATTGTACACATTCTTTGTTCAAGTGGTCTCCCGAT
GTCTATGGAAGATGCCAAAGACAAAGGCTTTGTGGTGGTGGAGAAGGAAGAACTGAACATGATTGAC
AACTTCTCAGTGAGCCAACAACCAAGAGCTGGGAGATCATCACTGTTGAAGAGGCAAAGCGCAGGAAG
AGCACATGCAGCTACTATGAAGACGAGGACGAAGAGGTGCTGCCTGCTCCTACGGCCACAGCGCGCTC
CTGGAGAATATGCACATCGAGCAGCTGGCCCGACGCCTTCTGCAAGGGTGAAGGGTATCCATGGAGA
CTGGCCTATAGCACGTTAGAGCACGGGACCAGCTTAAAGACGCTCTACCGAAAATCGGCATCACTAGAC
AGTCTGTCTATTGGTCAATCAAGATATGGATAATCAGATTTTTGGAGCATATGCAACTCATCCTTTC
AAGTTCAGTGACCACTATTATGGCACAGGCGAACTTTTCTCTACACATTAGCCCTCATTTTAAGGTC
TTTAAGTGGAGTGGAGAAAATTCATACTTTATCAATGGAGACATAAGTCTTTTAACTTGGTGGTGGGA
GGGGGACGATTTGGTTTATGGCTAGATGCTGATTTATACCAGGACGAAGCAACTCTGCAGCACTTTC
AATAATGATATTCTTTCCAAAAAGGAAGACTTCATAGTTCAGGATCTGGAGGTGTGGGCATTTGATGA
```

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT  
 TACAAGGATGACGACGATAAGGTTAAACGGCCGGC

Restriction Sites: SgfI-MluI

<b>Plasmid Map:</b>	□
<b>ACCN:</b>	NM_181782
<b>Insert Size:</b>	2829 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<a href="#">NM_181782.4</a>
<b>RefSeq Size:</b>	6457 bp
<b>RefSeq ORF:</b>	2829 bp
<b>Locus ID:</b>	135112
<b>UniProt ID:</b>	<a href="#">Q8NI08</a>
<b>Cytogenetics:</b>	6q22.31-q22.32
<b>Protein Families:</b>	Druggable Genome
<b>MW:</b>	106.2 kDa

**Gene Summary:**

Enhances the transcriptional activities of several nuclear receptors. Involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.

[UniProtKB/Swiss-Prot Function]

Transcript Variant: This variant (1) and variants 3 and 4 encode the same longest isoform (1).