

Product datasheet for **SC317068**

GSG1L (NM_001109763) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GSG1L (NM_001109763) Human Untagged Clone
Tag:	Tag Free
Symbol:	GSG1L
Synonyms:	PRO19651
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001109763, the custom clone sequence may differ by one or more nucleotides ATGAAGACTAGCCGCGCGGCCGAGCGCTCCTGGCCGTGGCCCTGAACCTGCTGGCGCTG CTGTTTCGCCACCACCGCTTCTCACCACGCACTGGTGCCAGGGCACGCAGCGGGTCCCC AAGCCGGGCTGCGGCCAGGGCGGGCGGCCAACTGCCCAACTCGGGCGCCAACGCCACG GCCAACGGCACCCGCCCGCCCGCCGCCCGCCGCCACCGCCTCGGGGAACGGC CCCCCTGGCGGCGCGCTCTACAGCTGGGAGACCGCGACGACCGCTTCTCTTCAGGAAT TTCCACACCGGCATCTGGTACTCGTGCAGGAGGAGCTCAGCGGGCTTGGTGAAAAATGT CGCAGCTTCATTGACCTGGCCCCGCATCGGAGAAAGGGTCTGTGGCTGTCTGTGGTC TCCGAGGTGCTGTACATTCTGTGCTGGTGGTTGGCTTACGCCTCATGTGTCTCGAGCTC TTCCACTCCAGCAATGTATCGACGGGCTCAAGCTCAATGCCTTCGCGGGTGTCTTCACG GTGCTCTCAGGCCTCCTGGGAATGGTCGCCACATGATGTACACGCAGGTGTTCCAGGTC ACCGTGAGCCTCGGTCTGAGGACTGGAGACCCATTCTGGGACTACGGGTGGTCTTCTC TGCTGGCGTGGGGTCTTTACCTGCTGCATGGCAGCCTCTGTACCACGCTCAACTCC TACACCAAGACGGTCATTGAGTCCGGCACAAGCGCAAGTCTTTGAGCAGGGCTACCGG GAAGAGCCGACCTTCATAGACCTGAGGCCATCAAGTACTTCCGGGAGAGGATGGAGAAG AGGGACGGGAGCGAGGAGGACTTTCACTTAGACTGCCGCCACGAGATACCTTGCCCGA CACCAGCCACACATGGCGGATTCTGGCCCCGGAGCTCCGCACAGGAAGCACCAGAGCTG AACCGACAGTGCTGGGTCTTGGGCACTGGGTG
Restriction Sites:	Please inquire
ACCN:	NM_001109763
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001109763.1, NP_001103233.1</u>
RefSeq Size:	4922 bp
RefSeq ORF:	996 bp
Locus ID:	146395
UniProt ID:	<u>Q6UXU4</u>
Cytogenetics:	16p12.1
Protein Families:	Transmembrane
Gene Summary:	<p>As a component of the inner core of AMPAR complex, modifies AMPA receptor (AMPA) gating.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (2) lacks several 5' coding exons and has an alternate 5' UTR, compared to variant 1. This results in a shorter protein (isoform 2) that lacks the GSG-1 conserved domain.</p>