

Product datasheet for SC317092

TRIM49D2 (NM_001105522) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TRIM49D2 (NM_001105522) Human Untagged Clone
Tag:	Tag Free
Symbol:	TRIM49D2
Synonyms:	TRIM48L1; TRIM49D1; TRIM49D2P; TRIM49L; TRIM49L1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317092 representing NM_001105522. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
 ATGAATTCTGGAATCTCGCAAGTCTTCCAGAGGGAACACCTGCCCATCTGCCTGAACTACTTCATA
 GACCCGGTCACCATAGACTGTGGGCACAGCTTTTGCAGGCCCTGTTTCTACCTCACTGGCAAGACATC
 CCAATTCTTACTCAGTGTCTTGAATGCCTAAAGACAACACAGCAGAGAAACCTCAAACTAACATTTCGA
 TTGAAGAAGATGGCTTCCCGTGCCAGAAAAGCCAGTCTCTGGCTATTCTGAGCTCTGAGGAGCAAATG
 TGTGGCACTCACAGGGAGACAAAGAAGATATTCTGTGAAGTGGACAGGAGCCTGCTCTGTTTGCTGTGC
 TCCAGCTCTCTGGAGACCGGTATCACAGACACTGTCCCGCTGAGTGGGCTGCTGAGGAACACCGGGAG
 AAGCTTTTAAAGAAAATGCAGTCTTTATGGGAAAAAGTTTGTGAAAATCAGAGAAACCTGAACGTGGAA
 ACCACCAGAATCAGCCACTGGAAGGATTATGTGAATGTAAAGCTAGAAGCTATTAGAGCTGAGTATCAG
 AAGATGCCTGCATTTTCATCATGAAGAAGAAAAACATAATTTGGAGATGCTGAAAAAGAAGGGGAAAGAT
 ATTTTTCATCGACTTCATTTAAGTAAAGCCAAAATGGCTCACAGGAGGGAGATTTTAAGAGGAACGTAT
 GCGGAGCTGATGAAAATGTGCCATAAACAGATGTGGAGCTACTTCAGGCTTTTGGAGACATATTACAC
 AGGAGTGAGTCCGTGCTGCTGCACATGCCCCAGCCTCTGAATCTAGAGCTCAGGGCAGGGCCCATCACT
 GGACTGAGGGACAGGCTCAACCAATTCGAGTGGATATTACTCTGCCTCATAATGAAGCCAACAGTCAT
 ATCTTCCGACGTGGAGATTTGAGAAGCATTGTATTGGATGTGACCGTCAAAATGCGCCCATATCACT
 GCAACACCTACAAGTTTCTTGCAATGGGTGCTCAGACTTTCACCTCTGGCAAAATATTACTGGGAGGTC
 CATGTGGGGACTCTTGAATTGGGCCTTTGGTGTCTGTAATAAGTATTGAAAGGGACGAATCAGAAT
 GGCAATATACATGGAGAGGAGGGACTCTTGTAGTCTGGGTGTGTTAAGAATGACATTCACTGCAGCCTC
 TTTACCACCTCCCACTTACACTGCAGTATGTCCAAGACCTACCAACCATGTAGGATTATTCCTGGAT
 TGTGAAGCTAGAACTGTGAGCTTCGTTGATGTTAATCAAAGCTCCCTATACACACCATCCCTAATTGC
 TCCTTCTACCTCCTCTCAGGCCTATCTTTGCTGTGTTACCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001105522
Insert Size:	1359 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105522.1</u>
RefSeq Size:	1651 bp
RefSeq ORF:	1359 bp
Locus ID:	729384
UniProt ID:	<u>C9J1S8</u>
Cytogenetics:	11q14.3
MW:	52.5 kDa