

Product datasheet for **SC317093**

KLHL33 (NM_001109997) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	KLHL33 (NM_001109997) Human Untagged Clone
Tag:	Tag Free
Symbol:	KLHL33
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317093 representing NM_001109997.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCACGATCGCC
ATGCTCCTGAGCGGGATGAGGGAATCCCAGGGCACAGAGGTATCTCTGCGGACGATCTCCACCCAGGAC
CTGCGACTCCTCGTCTCTTTTGTCTACTCCGAGTTGTGCGGGCAAGGTGGCCAGGGCTACTGAGAGCT
GCCCAGGCTGCTCTGCAGTACCAGAGCTCTTCTGCTTGGATTTGTGTCAGAAAGGCTTGGCACGGGGC
CTCAGCCCTGCCCGTTGCTGCGCCTGTTCCCATGGCGGAAGCCCTGGGTTGGAGAGGCTCTGGAGC
AAAGCCCGTCACTACCTCCTCACCCACCTGCCTGCTGTAGCCTTGTGCTGCTTTCCCTTCTTTACCA
GCTGCCTGCTTGGCTGAGCTCCTGGATAGTGATGAGCTCCATGTGCAGGAGGAGTTTGGGCCTTTGTG
GCTGCACGGTGTGGCTGGCTGCCAACCCGAGACCCAGGAGTCAGAGGCCAAGGCCCTGCTGCGATGT
GTCCGCTTTGGCCGCATGTCCACAGGGAGTTGCGGAGGGTGGGGCAGCCGGGCTACTTCCACCCCTG
ACCCAGATCTGTTGCACAGCTGATGGTAGAGGCTGATGTTCCAGGCCAAGAGAGACGGAGGGAGCCT
GACCGGGCACTGGTAGTGATTGGCGGGGATGGGCTCAGACCAGACATGGCCCTAAGACAACCATCCCGA
GCAAGTGTGGTGGGCCCGGCCCTTCCGCTGTGGCGTGGGACTGGTACGAACGTGTTAGTGGGGGCAGCTG
CCTGCCCTGCCTGCCCGGACGCTTCCGGCATGGGGCTGCAAGCCTGGCAGGAAGTGAACCTCTATGTG
TGTGGGGGACAAGATTTCTACAGTCACTCCAACACCCTGGCTTCAACTCTCAGGTGGGAGCCAGTCAA
GAGGACTGGGAGGAGATGGCTCCTTTGTCCAGGCTCGAAGCCTTTTCTCCTTGGTGGCACTGGATGGA
AACTTTATGCCCTGGGTGGAAGACACAATGATGTTGCCCTGGACTCTGTGGAGACCTACAACCCCTGAG
CTCAATGTCTGGAGGCCAGCACCTGCACCTCCAGCACCATGTTTGGCCACGCAGCTGCGATTTTGGAG
GGCCAGTTGTACGTGAGCGGTGGCTGTGGTGGGACTGGCCAATACCTGGCCTCACTGATGCACTATGAC
CCCAAACCTTGAGAAGCCAGGGACGTTTCTGAGCCCTATGGGGTACCTCGGGCTGGCCATGTCATGGCT
GCATTGGGTGGGAGGTTGTATGTGGCAGGTGGGCTGGGTGAGACTGAGGACCTGCTAAGCTTTGAGGCC
TATGAACCTAAGGACTGATAGCTGGACTCACCTGGCACCCTACCTCCCCCATGTGGGGCTGCAAGT
GCTGTGCTGCAGGGGAGCTACTGGTGTCTGGGGGCTACAGTACCCTACTTATGCCCTCTCTCACCTT
ATCCATGCCTACTGTCTGGCCTGGGCCGATGGCTCTGCCTGGAACTCTGCCAAGGCCTCGGGCTGAG
ATGCCTGCCTGCATCCTGACACTGCCCGCTGTGCAGCACATAGCTTTGGTTCCACCCACACCAAAAC
AAACCTGCTGGGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001109997

Insert Size: 1602 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: NM_001109997.2

RefSeq Size: 1863 bp

RefSeq ORF: 1602 bp

Locus ID: 123103

UniProt ID: A6NCF5

Cytogenetics: 14q11.2

MW: 57.8 kDa