

## Product datasheet for SC317099

### PRAMEF17 (NM\_001099851) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PRAMEF17 (NM_001099851) Human Untagged Clone
Tag:	Tag Free
Symbol:	PRAMEF17
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317099 representing NM_001099851. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAGCCTCCAGTCCCCATCCAGACTCCTGGAGCTGGCAGGCCAGAGCCTGCTGAGGAACCAAGTTCTTG
ACCATCTTCATCCTGGACGAGCTGCCAGGGAGGTCTTCCCTCTGATGTTTCATGGAGGCCTCCAGCATG
AGACATTTTGGAGCCCTGAAGCTGATGGTGCAGGCCTGGCCCTTCCCTCCGCCTCCCTCTGGGATCCCTG
ATGAAGACACCTCATCTGGAGACCTTGAAGCTGTGCTGAAGGGACTTGATACACTGTGGCCAGAAAG
CTTCGCCCCAGGAGGTGGAACCTCAAGTGTGGATTTGCGGGATGTTGATGGGAATTTCTGGACTATA
TGGTCTGGAGCCAGGCCCTCTCCTGCTCCCGAGGCCATGAGTAAGAGGCAGACAGTGGAGACTAT
CCAAGGACGGGAGAGACACCAGCCCTTGAAGGTGTTTCATAGACCTCTGCCAAAAGGAAAGTACACTGGAT
GAATGCCTGAGCTACCTCTGCAGGTGGATCCACTACAGAAGAGGTCTAGTGCACCTGTGTTGTAATAAG
GTGCAGAATTAATCAATGCCCACTTCAAGTTTCAGAAATTTATTGAAAAGGGTATACCCAGACAGTATC
CAGGAGTTGGAATTAAGAGAAAGTGTCTCTGAATAAAACAGGAAAGTTTGCCCTTACTTGAGCCAG
ATGAGCAATCTTCGAAACTCTTTTAGCCTTCGGTTATGACGATGAGTTATATGTAAGCGGCCAACAG
CAGTTCGTTCTGACTTGGACTGTCCATTCCTCTGCCTGACTACCCTCAGATGCTTTATATAAGAAAG
ATCAGTAATATCAAAGAGCACCTGGAGCACCTGCTCAGGTGCCTCAAGAACCCTTGGGAACCTTTATA
TTCTGTATGCTTACCTAGCTGATCAGGACATGGAGTGTCTGTCTCAGTACCCAAGCCTCAGTCAGCTA
AAGGAGCTGCATCTGATTCATATCCTAATGTGGACTACCAATCTTGAGCCCTTGGAGCTCTGCTAGAG
AAAGTTGCTGCTACTCTCGAGATCCTCACGTTAAAGGACTGTGAGATCCAGGACTCCCAGCTCAGGGTC
CTCCTGCCTGCCCTGAGCCGCTGCTCCCAGCTCACCACCTTCTACTTTGCGGAAATGAGACCTCCACG
AATGCTCTGAAAGACCTGCTGTGTACACAGGTTGGGCTGAGCAAGTTAGGTCTGGAGTTGTATCCTGCC
CCTCTGGAGTGTCTTGACAACAGGGTGTCAATTGGGAGATCCTCGCCCAATTGCGGCTGAGCTG
ATGTGTACTCAGGGAAGTCAGGCAGCCCAAGAGGATCTTTTTTGGTCCCATCCCCTGCCCTTCTGT
GGCTCATGGCCATCTGAGAAAGTGGACTTCCATCTTTGCTCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```



[View online >](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001099851
<b>Insert Size:</b>	1425 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001099851.3</a></u>
<b>RefSeq Size:</b>	1578 bp
<b>RefSeq ORF:</b>	1425 bp
<b>Locus ID:</b>	391004
<b>UniProt ID:</b>	<u><a href="#">Q5VTA0</a></u>
<b>Cytogenetics:</b>	1p36.21
<b>MW:</b>	54.8 kDa