

Product datasheet for **SC317103**

ZNF83 (NM_001105554) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF83 (NM_001105554) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF83
Synonyms:	HPF1; ZNF816B
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001105554, the custom clone sequence may differ by one or more nucleotides ATGCATGGGAGAAAGGATGATGCACAAAAGCAGCCTGTTAAAAATCAGCTTGGATTAAC CCGCAGTCACATCTACCAGAAGTGCAGCTATTTCAAGCTGAAGGGAAAATATATAAATAT GATCACATGGAAAAATCTGTCAACAGTAGTTCCTTAGTTTCCCCACCCCAACGTATTTCT TCTACTGTCAAAACCCACATTTCTCATAACATGAATGTAATTTTGTGGATTTCATTATTC ACACAAAAAGAGAAAGCAAATATTGGGACAGAACTACAATGTAGTGAGCGTGGCAAG GCCTTTCATCAAGGCTTACATTTACTATACATCAAATAATCCATACTAAAGAGACGCAA TTTAAATGTGATATATGTGGCAAGATCTTCAATAAAAAATCAAACCTTGCAAGTCATCAA AGAATTCATACTGGAGAGAAGCCATATAAATGTAATGAATGTGGCAAGGCTTCCATAAT ATGTCACACCTTGACAGCATCGCAGGATTCATACTGGAGAGAAAACCATATAAATGTAAT GAATGTGGCAAGGCTTTAATCAAATTTACACCTTGACAACATCAAAGAATTCATACC GGAGAGAAACCTTATAAATGTAATGAATGTGGAAAGGCTTCCATCAAATTTACACCTT GCACAACATCGGACAATTCATACTGGAGAAAAACCTTACGAATGTAACAAATGTGGCAAG GTGTTCAAGTCGCAATTCCTACCTTGTACAACATCTGATCATTACTGGAGAGAAACCT TACAGATGTAATGTATGTGGAAAGGCTTCAGTCACAAGTCATCCCTAGTAAATCACTGG AGAATTCATACTGGAGAGAAACCTTACAATGTAATGAGTGTGGCAAGGCTTTCAGTCAC AAGTCATCCCTAGTAAATCACTGGAGAAATCCACACTGGAGAGAAACCTTACAATGTAAT GAATGTGGCAAGGCTTTCAGTCGCAATTCATACCTTGCCCAACATCTGATAATTCATGCC GGTGGAGAAACCTTATAAGTGTGATGAATGTGACAAAGCATTTCAGTCAAAATTCACATCTT GTACAACATCACAGAATCCATACTGGAGAGAAACCTTACAATGTGATGAATGTGGCAAA GTCTTCAGTCAAAATTCATACCTTGCATATCATTGGAGAATTCATACTGGAGAAAAAGCT TATAAATGTAATGAATGTGGGAAGGCTTCGGTCTAAACTCATCCCTAGCACATCATCGG AAAATTCACACTGGAGAGAAACCTTCAAATGTAATGAATGTGGCAAAGCTTTTAGTATG CGTTCAAGCCTCACTAATCATCATGCGATCCACACTGGAGAGAAACATTTCAAATGTAAT GAATGTGGCAAACCTTCCGCGACAATTCATATCTGTACGTCATCAGAGATTTTCATGCC GGAAAGAAATCTAACACATGTAAT
Restriction Sites:	Please inquire
ACCN:	NM_001105554



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105554.1</u> , <u>NP_001099024.1</u>
RefSeq Size:	2538 bp
RefSeq ORF:	1467 bp
Locus ID:	55769
Cytogenetics:	19q13.41
Protein Families:	Transcription Factors
Gene Summary:	<p>May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]</p> <p>Transcript Variant: This variant (7) lacks three consecutive exons in the 5' UTR and lacks an in-frame segment in the CDS, compared to variant 1. The resulting isoform (b) lacks an internal segment, compared to isoform a.</p>