

Product datasheet for **SC317114**

USP17 (NM_001105662) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	USP17 (NM_001105662) Human Untagged Clone
Tag:	Tag Free
Symbol:	USP17
Synonyms:	MGC119330; MGC119331; MGC119333; RS447; USP17A; USP17H; USP17I; USP17J; USP17K; USP17L; USP17M
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317114 representing NM_001105662.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGGACGACTCACTCTACTTGGGAGGTGAGTGGCAGTTCAACCCTTCCAAAACCTCACATCTTCT
CGGCCCGATGCAGCTTTTGTGAAATCCAGCGGACTTCTCTCCCTGAGAAGTCACCCTCATGTGAG
ACCCGTGTCGACCTCTGTGACGATTTGGCTCCTGTGGCAAGACAGCTTGCTCCCAGGGAGAAGCTTCT
CTGAGTAGCAGGAGACCTGCTGCCGTGGGGCTGGGCTCCAGAATATGGGAAATACCTGCTACGTGAAC
GCTTCCTTGGAGTGGCTGACATACACACCGCCCTTGCCAACTACATGCTGTCCCGGGAGCACTCTCAA
ACGTGTCACTCGTCAAGGGCTGCATGCTGTACTATGCAAGCTCACATCACACGGGCCCTCCACAAT
CCTGGGCACGTATCCAGCCCTCACAGGCATTGGCTGCTGGCTTCCATAGAGGCAAGCAGGAAGATGCC
CATGAATTTCTCATGTTCACTGTGGATGCCATGAAAAAGGCATGCCTTCCCGGGCACAAGCAGGTAGAT
CATCACTTAAGGACACCACCTCATCCACCAATATTTGGAGGCTACTGGAGATCTCAAATCAAGTGT
CTCCACTGCCACGGCATTTCAGACACTTTTGACCTTACCTGGACATCGCCCTGGATATCCAGGCAGCT
CAGAGTGTCCAGCAAGCTTTGGAACAGTTGGTGAAGCCGGAAGAACTCAATGGAGAGAATGCCTATCAT
TGTGGTGTGTGCTCCAGAGGGCGCCGGCTCCAAGACGTTAACTTTACACACCTCTGCCAAGGTCCCTC
ATCCTTGTATTGAAGAGATTCTCCGATGTCACAGGCAACAAGATTGCCAAGAATGTGCAATATCCTGAG
TGCTTGCATGCAGCCATACATGTCTCAGCAGAACACAGGACCTTGTCTATGTCTCTATGTCTGTG
CTGGTCCACGCTGGGTGGAGTTGTCAACGGACATTACTTCTTATGTCAAAGCTCAAGAAGGCCAG
TGGTATAAAATGGATGATGCCGAGGTCACCGCTCTAGCATCACTTCTGTCTGAGTCAACAGGCCCTAC
GTCCCTTTTACATCCAGAAGAGTGAATGGGAAAGACACAGTGAAGTGTGTCAAGAGGCCAGGAACCA
AGAGCCCTTGGCGAGAAGACACAGACAGCGGAGCAACGCAAGGAGAGCTCAAGAGAGACCACCCCTGC
CTCCAGGCCCCGAGTTGGACGAGCACTTGGTGGAAAGAGCCACTCAGGAAAGCACCTTAGACCCTGG
AAATTCCTTCAAGAGCAAAACAAAACGAAGCCTGAGTTCAACGTGAGAAAAGTCAAGGTACCCTGCCT
CCCAGCTACTTGTGATTCATCAATCAAATACAAGTGTGGGATGAAGAACCATCATCTGAACAGCAA
AGCTCCCTGTAAACCTCTCTTCGACGACCCCGACACATCAGGAGTCCATGAACACTGGCACACTCGCT
TCCCTGCGAGGGAGGGCCAGGAGATCCAAGGGGAAGAACAACACAGCAAGAGGGCTCTGCTGTGTGC
CAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001105662

Insert Size: 1593 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001105662.1](#)

RefSeq Size: 1593 bp

RefSeq ORF: 1593 bp

Locus ID: 391627

Cytogenetics: 4p16.1

MW: 59.8 kDa