

Product datasheet for **SC317122**

LINGO3 (NM_001101391) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LINGO3 (NM_001101391) Human Untagged Clone
Tag:	Tag Free
Symbol:	LINGO3
Synonyms:	LERN2; LRRN6B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317122 representing NM_001101391.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGC
ATGACCTGCTGGCTGTGCGTCTGAGCCTGCCCTGCTCCTGCTGCCCGCGCGCCCGCCCGGCTGGA
GGCTGCCCGGCCGCTGCGAGTGCACCGTGCAGACCCGCGCGTGGCTGCACGCGCCGCGCTGACC
GCGTGCCCGACGCGATCCCGGCCGAGACCCGCTGCTGGAGCTCAGCCGCAACCGCATCCGCTGCCTG
AACCCGGGCGACCTGGCCGCGCTGCCGCGCTGGAGGAGCTGGACCTGAGCGAGAACGCCATCGCGCAC
GTGGAGCCCGGCCCTTCGCAACCTGCCGCGCTGCGCGTCTCGCTCTCCGTGGCAACAGCTGAAG
CTCATCCCGCCGGGTCTTCACGCGCTGGACAACTCACGCTGCTGGACCTGAGCGAGAACAGCTG
GTAATCTGCTGGACTACACTTTCAGGACCTGCACAGCTGCGCCGGCTGGAAGTGGGCGACAACGAC
CTGGTATTCGCTCGCGCCGCGCTTCGCGGGCTGCTGGCCCTGGAGGAGCTGACCTGGAGCGCTGC
AACCTCACGGCTCTGTCGGGGAGTCTGCGGCCATCTGCGCAGCTGGCGCCCTGCGGCTGCGCCAC
CTGGCCATCGCTCCCTGGAGGACAGAACTTCCGAGGCTGCCCGGCTGCTGCACCTGGAGATTGAC
AACTGGCCGCTGCTGGAGGAGTGGCGGGGCGAGCCTGCGGGGCTGAACCTGACCTCGCTGTCGGTC
ACCCACACCAACATCACCGCCGTGCCGGCGCCGCGCTGCGGCACAGGCGCACCTCACCTGCCTCAAT
CTGTGCGACAACCCATCAGCACGGTGCCGCGGGGTCTTCCGGGACCTGGTCCGCTGCGCGAGCTG
CACCTGGCCGGGGCCCTGCTGGCTGTGGTGGAGCCGAGGCTTCTGGGCTGCGCCAGATCCGCTG
CTCAACCTCTCCAACAACCTGCTCTCACGTGGAGGAGAGCACCTTCCACTCGGTGAACACGCTAGAG
ACGCTGCGCGTGGACGGGAACCGCTGGCTGCGACTGTCGCTGCTGTGGATCGTGCAGCTCGCAAG
ACCCTCAACTTCGACGGGCGGCTGCCGGCTGCGCCACCCGGCCGAGGTGCGCGGCGACGCGCTGCGA
AACCTGCCGACTCCGTGCTGTTGAGTACTTCTGTGCGCAAAACCAAGATCCGGGAGCGGCGCTG
CAGCGCGTACCGCCACCGGGGCGAAGACGTCCGCTTCTCTGCGCGCCGAGGGCGAGCCGGCGCCC
ACCGTGGCTGGGTGACCCCCAGCACCGGCGGTGACGGCCACAGCGGGGCGGGCGCGCTGCTC
CCCGGGGGGACGCTGGAGATCCAGGACGCGCGCCGAGGACAGCGGCACCTACAGTGCCTGGCAGC
AACGCGGGCGGCAACGACACCTACTTCGCCACGCTGACCGTGCGCCCGAGCCGGCGCCAAACGGACC
CCGGGCGAGGCCACAACGAGACGCTGGCGGCCCTGCGCGCGCCGCTCGACCTACCAACCATCTGGTG
TCCACCGCATGGGTGCATCACCTTCTGGGCGTGGTCTTCTGCTTCTGCTGCTGTTCTGTGTGG
AGCCGCGCGCGGGCAGCACAAAAACAACTTCTCGGTGGAGTACTCTTCCGAAGGTGGATGGGCCG
GCCGCGCGGGCGGCGAGGAGCGCGCAAGTTCAACATGAAGATGATTGA
ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_001101391

Insert Size: 1779 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001101391.1](#)

RefSeq Size: 2130 bp

RefSeq ORF: 1779 bp

Locus ID: 645191

UniProt ID: [P0C6S8](#)

Cytogenetics: 19p13.3

Protein Families: Transmembrane

MW: 64.9 kDa