

Product datasheet for **SC317123**

LRRIQ3 (NM_001105659) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRIQ3 (NM_001105659) Human Untagged Clone
Tag:	Tag Free
Symbol:	LRRIQ3
Synonyms:	LRRC44
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317123 representing NM_001105659.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTTTCATGGAACAGTCACAGAAGAGCTAACCGTCAATGAAGAATGGAGTCACTATAATGAAAACATA
AGAGAAGTCAAAAAGATTTTGTGTTTGTGAAGTCAATGGCCTTCATTTAAAGTCTATGGAAAATTTG
CAGTCTTGCACTCTCTTAGAGTATGCATCTTCTCAAACAATTTTATTACAGATATTCACCACTCAA
AGTTGTATAAAATTAATCAAACCTGATCTCCATGGAAATCAGATAAAGAGTCTACCAAATACCAAATTT
TGAATGGATTGAAGAACCTAAACTACTCTATCTTCATGACAATGGGTTTGCAAAGTTAAAGAATATA
TGTGTATTATCTGCCTGTCCAACCTCATTGCCCTCACTATGTTTGATTGTCCAGTAAGCCTTAAAAAA
GGATATAGACATGTTCTTGTTAACAGTATATGGCCTCTCAAAGCGCTGGATCATCATGTGATTTCTGAT
GAAGAAAATAATTCAGAACTGGCATCTTCTGAAAGATTCAAAGCATGTAACCATCGACTTTTCTTTAAT
TTCTGCCAGCTTTGAGAAAGGGAACAACCTATGAAGAGGAAATTAATAATTAACATATTACTTCA
AAAATTAATGCAATTCTGGCTCATAATTCACAGTTTTGATTGTTCAAAGATGGATACGTGGTTTCTTA
GTTAGAAAAAATTTGAGCCCTGTGTTTTCCACAAAAAAAACAGCAGGAAAAAATATTAGAGGATAT
GAAGCAAATGGATTTACATAACCAAAGGGTATGAAGATAAGCTCCTTAAGGATCTCTTTTTCAAACCT
GAAACTAATATAAAAGGAAAGCTTGCAATTTGGAAACATAATATATATTCTGTGTTTAAAAAAT
TCCAGTGAACACAGAAAACATGTGTCACTATTTTATGTGAATTAACCAAAGATCTAGGCATGAAA
TCAAAAACATCAAGACATCTCATTCAAAAAGGTCAGGAGTCTGAAGATGAAATTTGGATGAAAAATG
GATACCAGTTTTAGGATATCAGTTTTCAAACCTACCATATACACTTCAGGTTTCAATGAAGAATAATGCA
GTATTGAGAGAGAAAAACAACATTTTTTCTGCATATCCTCAGCCAATCTATACTACTCATCCAAAG
CCAATCATTA AAAAGACATACGATTGGAGCGGAGTATGAAAGATTTTTTGCACCACAAAGAGCTGGT
ATGAAACTCCGAACATTTAGTGATATTGACAAATATTACACAGAACAAAAGAGCAGGAATACCATAAA
GAAAAAGTAAGAGTTGTAGCCATGGCACAAAGTTGCTCGAGAAAAGAGTTAGAGTAGCTGTTAATGAACAT
TTGAATCAGAAAAATATGCTACACAAAACTAATTGAAGAAAATAAAGAGACAATTCAGAACAGTTTA
CGACAAGTTTGGCAAAACAGATTCAACTACCTTAAAAAGCTAGAGAGAGAAAAGCTCTGTTTCTAAAA
GAAAAATCCCAAAGGCTTCAGAGCGTTTATTAGTTCAAACCTAAATAATGAGCGCACTTTTTGACC
AGAGGACTACTTAAAATGACAGACTGGAGAAGAATGAGGCTGTCCTAAAAGAGAAAAGCCTGATTGTT
AAGCAAAAATAAAGCAGAAAAATATAGAAAGAAATTTACTTAAAGAAAATAAAGAAAGTTAGATCTCAA
GAAATATATAAAAGACATTGTGAAGAAAATTTGTTATGGATATGATTGCCTTTGAAAAAGCCTGTGAA
AGACTTCAAGATGCTAAAACAAAAGTAGCAATTGTGAAAACAAATTTAGACTTTAAAGTTCCAATGGA
CTGATAAAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: Sgfl-Mlul

Plasmid Map: □

ACCN: NM_001105659

Insert Size: 1875 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105659.1</u>
RefSeq Size:	2865 bp
RefSeq ORF:	1875 bp
Locus ID:	127255
UniProt ID:	<u>A6PVS8</u>
Cytogenetics:	1p31.1
MW:	73.7 kDa