

Product datasheet for **SC317127**

WDR43 (NM_015131) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	WDR43 (NM_015131) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR43
Synonyms:	NET12; UTP5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317127 representing NM_015131.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGCACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGCGGCGGGCGGGCGGGTAGCTGCGACCCCTGGCCCTGCTGGGGTCCCTTGCCTTCTCCCCG
CACAGCCAGGCCTACTTCGCTTTGGCCTCTACCGACGGTCACTTACGAGTATGGGAGACGGCCAACAAC
CGGCTGCACCAAGGAGTACGTGCCTTCCGCGCACCTCAGTGGTACCTGCACCTGTCTGGCCTGGGCGCCA
GCGCGGCTGCAGGCCAAGGAAAGTCCCCAGAGGAAAAAAGGAAATCAGAAGCTGTAGGAATGAGTAAC
CAGACTGACTTATTGGCTCTTGGCACAGCAGTTGGTAGCATTTTATTATACAGCACAGTAAAGGAGAG
TTACACAGTAAATTAATAAGTGGTGGACATGACAACAGAGTCACTGCATACAGTGGCATCAAGACAGT
GGCTGTTTATATAGTTGTTTCAGATGATAAATATTGTGGAATGGAACGTACAGACATGCAAGTAAAG
TGCAAATGGAAAGGCGACAATAGCAGTGTCAAGTCCCTATGTATCAGCCAGATGGAAAGATGTTGCTT
TCAGCTGGTGAACAATCAAATATGGGTTTTGGAGACCAAGAAGTCTACAGGCATTTACAGGACAT
GCAACGCCAGTTTCGTCACTGATGTTCACTACCATCAGACCTCCTAATGAGAGCCAGCCCTTTGATGGA
ATTACAGTCTTTATTTCTTATCTGGAGCAGTACATGACCGGTTACTTAATGTCTGGCAGTCCGATCA
GAAAACAAAGAAAAGAGTGCAGTGTATTTACAGTTACCGATGAACCTGTCTATATTGACTTAACT
TTGTGAGAAAACAAAGAGAGCCTGTCAAGTTGGCTGTTGTTTGCAGAGATGGTCAAGTCCATCTTTT
GAACACATATTAAATGGGTACTGCAAAAAGCCTTTGACTTCAAATGCACAATTCAGATAGCAACACCT
GGGAAAGGCAAGAAGTCAACACCAAAACCCATCCCTATTCTAGCTGCTGGTTTTTGTCTAGACAAAATG
TCATTGTTGCTTGTATATGGCAGTTGGTTTCAGCCTACTATTGAGCGAGTGGCTTTAACTCCAGAGAA
CCTCATATGTGTTTAGTAAGAGATATTTCAAAGTCTGGGCCCCAAAGTAGAAACAGCTATAACAAAG
GTGAGGACACCAAGTGAATTTCTGAAGCAAAAGTTCTGGTGCCTGGGATTCTGGTGCATCATGCAGCT
ATCAAGCCCGCTCTCCACAAACCGAGCAAGTAGAGAGCAAGAGGAAGTCAGGGGGAATGAGGTTAGC
ATTGAAGAACGTCTGGGAGCAATGGATATAGACACACAAAAAAGGAAAGGAAAGACCTCCAGACGAAT
AGCTTTCCAGTCTTCTTACCCAGGGCTTAGAAAGTAACGATTTTGAAATGCTAAATAAAGTACTTCAA
ACTAGGAATGTAAACCTTATAAAGAAGACTGTATTAAGGATGCCCTGCATACTATTATCCGTTGTTA
CAAGAGCTTACAAAGAGGTTACAAGGACATCCTAATAGTGTGTGCTAATGGTTCAGTGGCTAAATGT
GTGTTAACAGTTCATGCATCATACTGTCCACGTTGCCTGACCTGGTACCCAGCTGGGACACTCTAC
CAGTTAATGGAAAGCAGAGTCAAACTTTTCAGAACTTTACACCTTCATGGAAAGCTTATTCTTCTA
ATTACACAAGTAACAGCATCAGAGAAGACAAAGGAGCAACTTCCCTGGACAGAAGGCAAGTTGGTG
TATGAAGAAGAGTCTTCTGAAGAGGAGTCTGATGATGAAATAGCAGATAAGGATTCTGAAGATAATTGG
GATGAAGATGAGGAGGAGAGTGAAGTGAAGGAGATGAGGACGTTGAAGAGGAAGATGAGGATGCCGAA
GGAAAAGATGAAGAAAATGGCGAGGACAGAGATACAGCAAGTGAAGGAAATTAATGGAGATTCTGAT
TTAGATCCTGAAAATGAAGTGAAGAAGATGA
ACGCGTACGCGGCCGCTCGAGCGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
  
```

Restriction Sites: SgfI-MluI

ACCN: NM_015131

Insert Size: 2034 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_015131.2</u>
RefSeq Size:	3542 bp
RefSeq ORF:	2034 bp
Locus ID:	23160
UniProt ID:	<u>Q15061</u>
Cytogenetics:	2p23.2
MW:	74.9 kDa
Gene Summary:	Ribosome biogenesis factor. Involved in nucleolar processing of pre-18S ribosomal RNA. Required for optimal pre-ribosomal RNA transcription by RNA polymerase I.[UniProtKB/Swiss-Prot Function]