

Product datasheet for **SC317130**

TLE3 (NM_001105192) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TLE3 (NM_001105192) Human Untagged Clone
Tag:	Tag Free
Symbol:	TLE3
Synonyms:	ESG; ESG3; GRG3; HsT18976
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

Fully Sequenced ORF: >SC317130 representing NM_001105192.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGGATCGCC
ATGTATCCGCAGGGCAGACATCCGGTCCCATCAACCCGGGCAGCCGGGATTTAAATTCACGGTGGCT
GAGTCTTTGTGACAGGATCAAAGACGAATTCAGTTCCTGCAAGCTCAGTATCACAGCCTCAAAGTGGAG
TACGACAAGCTGGCAAACGAGAAGACGGAGATGCAGCGCCATTATGTGATGTAATGAGATGTCCTAT
GGCTTGAACATTGAAATGCAACAAGCAGACAGAGATTGCGAAGAGACTGAACACAATTTTAGCACAGATC
ATGCCTTTCTGTCAACAAGACACCAGCAGCAGGTGGCGCAGGCAGTGGAGCGGCCAAGCAGGTCAAC
ATGACGGAGCTGAACGCCATCATCGGCAGCAGCAGCTCCAGGCGCAGCACCTCTCCATGCCACACAC
GGCCCCCGGTCCAGTTGCCACCCACCCGTCAGGTCTCCAGCCTCCAGGAATCCCCCAGTGACAGGG
AGCAGCTCCGGGCTGCTGGCACTGGCGCCCTGGGCAGCCAGGCCATCTGACGGTGAAGGATGAGAAG
AACCACATGAACTCGATCACAGAGAGAGAATCCAGTGCGAATAACTCTGTGTACCCCTCGAAAGC
CTCCGGGCCAGTGAGAAGCACCGGGGCTCTGCGGACTACAGCATGGAAGCCAAGAAGCGGAAGCGGAG
GAGAAGGACAGCTTGAGCCGATACGACAGTGATGGAGACAAGAGTGATGATCTGGTGGTGGATGTTCC
AATGAGGACCCCGCAACGCCCGGGTCAAGCCCGGCACACTCCCCTCCTGAAAATGGGCTGGACAAGGCC
CGTAGCCTGAAAAAGATGCCCCACCAGCCCTGCCTCGGTGGCCTCTTCCAGTAGCACACCTTCTCC
AAGACCAAGACCTTGGTCATAACGACAAATCCTCCACCCTGGGCTCAAGTCCAACACACCAACCCCA
AGGAACGACGCCCAACTCCAGGCACCAGCAGCAGCCAGGGCTCAGGTGATGCCGGTAAACCTCCG
GGCATGGACCCGATAGCCTCGGCTCTGCGCAGCCCATCTCCATCACCAGCTCCTATGCGGCGCCCTC
GCCATGATGAGCCACCATGAGATGAACGGCTCCCTCACCAGTCTGGCGCCTACGCCGCTCGCAATGGT
ATCCCACCCAGATGAGCGCCCGCGCTGCTGCAGCCGCTGCCTATGGCCGATCGCCAAATGGTGAAG
TTTGGAGCTGTTGGTTTTGACCCTCACCCCGATGCGGGCCACAGGCCCTCCCCTCAAGCCTGGCCTCC
ATTCTGGAGAAAACCAGCGTACTCATTCCATGTGAGTGTGATGGGCAGATGCAGCCGTGCCCTTC
CCCCACGACGCCCTGGCAGGCCCGGCATCCCAGGCACGCCCGGCAGATCAACACACTAGCCACGGG
GAGGTGGTGTGTGCCGTGACCATCAGCAACCCACGAGGCACGTCTACACAGGTGGCAAGGGCTGCGTG
AAGATCTGGGACATCAGCCAGCCAGGCAGCAAGAGCCCATCTCCAGCTGGACTGCCTGAACAGGGAC
AATTACATCCGCTCCTGCAAGCTGCTCCCTGATGGGCGCAGCTCATCGTGGGCGCGAGGCCAGCAG
CTCACCATCTGGACCTGGCCTCGCCACGCCCGCATCAAGGCCAGCTGACGTCTCGGCTCCCGCC
TGTTATGCCCTGGCCATTAGCCCTGACGCCAAAGTCTGCTTCTCCTGCTGCAGCGATGGGAACATTGCT
GTCTGGGACCTGCACAACCAGACCCTGGTCAGGCAGTTCAGGGCCACACAGATGGGGCCAGCTGCATA
GACATCTCCCATGATGGCACAAACTGTGGACAGGGGGCCTGGACAACACGGTGGCCTCCTGGGACCTG
CGGGAGGGCCGACAGCTACAGCAGCATGACTTCACTTCCCAGATCTTCTCGTGGGCTACTGCCCCACT
GGGGAGTGGCTGGCTGTGGGCAATGGAGAGCAGCAACGTGGAGGTGCTGCACCACCAAGCCTGACAAG
TACCAGCTGCACCTGCACGAGAGCTGCGTGTCTCCCTCAAGTTCGCCTACTGCGCAAGTGGTTCGTG
AGCACTGGGAAAGATAACCTTCTCAACGCCCTGGAGGACGCCTTATGGAGCCAGCATATCCAGTCTAAA
GAATCCTCGTCTGTCTTGAGTTGTGACATTTAGCGGATGACAAATACATTGTAACAGGCTCTGGTGAC
AAGAAGGCCACAGTTTATGAGGTCATCTACTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```

Restriction Sites: SgfI-MluI

Plasmid Map: □

ACCN: NM_001105192

Insert Size: 2310 bp

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001105192.2</u>
RefSeq Size:	5751 bp
RefSeq ORF:	2310 bp
Locus ID:	7090
UniProt ID:	<u>Q04726</u>
Cytogenetics:	15q23
Protein Families:	Transcription Factors
MW:	83.1 kDa
Gene Summary:	<p>This gene encodes a transcriptional co-repressor protein that belongs to the transducin-like enhancer family of proteins. The members of this family function in the Notch signaling pathway that regulates determination of cell fate during development. Expression of this gene has been associated with a favorable outcome to chemotherapy with taxanes for ovarian carcinoma. Alternate splicing results in multiple transcript variants. Additional alternatively spliced transcript variants of this gene have been described, but their full-length nature is not known. [provided by RefSeq, Sep 2013]</p> <p>Transcript Variant: This variant (2) uses an alternate in-frame splice site in the coding region compared to variant 1. The encoded isoform (b) is shorter than isoform a. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p>