

Product datasheet for SC317222

LYRM5 (ETFRF1) (NM_001001660) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LYRM5 (ETFRF1) (NM_001001660) Human Untagged Clone
Tag:	Tag Free
Symbol:	LYRM5
Synonyms:	LYRM5
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317222 representing NM_001001660. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAAAATGGCCAATTCTTTAAGAGGAGAAGTACTAAAACCTTTATAAAAACTGCTGTATCTTGGACGA
GACTATCCAAAAGGAGCAGACTATTTTAAAAAGCGTTTGAAGAACATTTTCCTTAAAAACAAGATGTG
AAGAAATCCAGAGAAGATCAAAGAAGCTATTGCACAGGGCGAATTTGTAATGAAAGAGCTAGAAGCTTTG
TACTTCCTTAGGAAATACAGAGCTATGAAACAACGCTATTATTCAGATACCAACAAAATAATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001001660
Insert Size:	273 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001001660.2</u>
RefSeq Size:	1149 bp
RefSeq ORF:	273 bp
Locus ID:	144363
UniProt ID:	<u>Q6IPR1</u>
Cytogenetics:	12p12.1
MW:	10.9 kDa
Gene Summary:	Acts as a regulator of the electron transfer flavoprotein by promoting the removal of flavin from the ETF holoenzyme (composed of ETFA and ETFB).[UniProtKB/Swiss-Prot Function]