

Product datasheet for **SC317231**

Trefoil Factor 3 (TFF3) (NM_003226) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Trefoil Factor 3 (TFF3) (NM_003226) Human Untagged Clone
Tag:	Tag Free
Symbol:	Trefoil Factor 3
Synonyms:	ITF; P1B; TFI
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_003226 edited AGATCTAACTGTCTGGGAGCTTGACAAAGGCATGCAGGAGAGAACAGGAGCAGCCACAGC CAGGAGGGAGAGCCTTCCCCAAGCAAACAATCCAGAGCAGCTGTGCAAACAACGGTGCAT AAATGAGGCCTCCTGGACCATGAAGCGAGTCTGAGCTGCGTCCCGGAGCCCACGGTGGT CATGGCTGCCAGAGCGCTCTGCATGCTGGGGCTGGTCCTGGCCTTGCTGTCTCCAGCTC TGCTGAGGAGTACGTGGGCTGTCTGCAAACAGTGTGCCGTGCCAGCCAAGGACAGGGT GGACTGCGGCTACCCCCATGTCACCCCCAAGGAGTGCAACAACCGGGGCTGCTGTTTGA CTCCAGGATCCCTGGAGTGCCTTGGTGTTCAGCCCTGCAGGAAGCAGAATGCACCTT CTGAGGCACCTCCAGCTGCCCGGGCGGGGATGCGAGGCTCGGAGCACCTTGCCCGG CTGTGATTGCTGCCAGGCACTGTTCATCTCAGCTTTTCTGTCCCTCTAGA
Restriction Sites:	Please inquire
ACCN:	NM_003226
Insert Size:	500 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_003226.2, NP_003217.2</u>
RefSeq Size:	685 bp
RefSeq ORF:	393 bp
Locus ID:	7033
UniProt ID:	<u>Q07654</u>
Cytogenetics:	21q22.3
Domains:	PD
Protein Families:	Secreted Protein
Gene Summary:	Members of the trefoil family are characterized by having at least one copy of the trefoil motif, a 40-amino acid domain that contains three conserved disulfides. They are stable secretory proteins expressed in gastrointestinal mucosa. Their functions are not defined, but they may protect the mucosa from insults, stabilize the mucus layer and affect healing of the epithelium. This gene is expressed in goblet cells of the intestines and colon. This gene and two other related trefoil family member genes are found in a cluster on chromosome 21. [provided by RefSeq, Jul 2008]