

Product datasheet for SC317251

C19orf12 (NM_001031726) Human Untagged Clone

Product data:

| | |
|---------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | C19orf12 (NM_001031726) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | C19orf12 |
| Synonyms: | MPAN; NBIA3; NBIA4; SPG43 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC317251 representing NM_001031726. Blue=Insert sequence Red=Cloning site Green=Tag(s) |

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGGAGAGGCTGAAGTCACACAAGCCCGCCACGATGACTATCATGGTGGAGGACATCATGAAGCTGCTG
TGCTCCCTTTCTGGGAGAGGAAGATGAAGCGGCTGTCAAGCACTCTGGGAAGGGTGCCCTGGTCACA
GGGGCCATGGCCTTCGTCGGGGTTTGGTGGGCGGCCACCGGGACTCGCCGTTGGGGGGCTGTCCGG
GGGCTGTAGGTGCTGGATGACAAGTGGACAGTTTAAGCCGTTCTCAGATCCTAATGGAGCTGCC
CCTGCCGAGCAACAGAGGCTCTTAACGAAGCCGAGCCATCATCAGGCACCTGGAGTGGACGGACGCC
GTGCAGCTGACCGCGCTGGTCATGGGCAGCGAGGCCCTGCAGCAGCAGCTGCTGGCCATGCTGGTGAAC
TACGTCACCAAGGAGCTGCGGGCCGAGATCCAGTATGATGACTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-Mlul

Plasmid Map:

ACCN: NM_001031726

Insert Size: 459 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

| | |
|-------------------------------|---|
| OTI Annotation: | This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001031726.3 |
| RefSeq Size: | 4412 bp |
| RefSeq ORF: | 459 bp |
| Locus ID: | 83636 |
| UniProt ID: | Q9NSK7 |
| Cytogenetics: | 19q12 |
| Protein Families: | Transmembrane |
| MW: | 16.3 kDa |
| Gene Summary: | <p>This gene encodes a small transmembrane protein. Mutations in this gene are a cause of neurodegeneration with brain iron accumulation-4 (NBIA4), but the specific function of the encoded protein is unknown. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011]</p> <p>Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.</p> |