

Product datasheet for SC317255

YPEL3 (NM_031477) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	YPEL3 (NM_031477) Human Untagged Clone
Tag:	Tag Free
Symbol:	YPEL3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317255 representing NM_031477. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGTGTGTGGCCAGGTCCTGACAGCCACCTACTCCCTCCCGGCAGGCACTGGGCTCCCTCTGCTCC
 CCGTGGGCGCTCCCGCGTGGGGCCACTGCCCGGCCCGCCCATGGTGCGGATTTCAAAGCCCAAG
 ACGTTTCAGGCCTACTTGGATGATTGTACCGGAGGTATAGCTGTGCCCACTGCCGCGCTCACCTGGCC
 AACCACGACGACCTCATCTCCAAGTCCTCCAGGGCAGTCAGGGGCGTGCCTACCTCTTCAACTCAGTG
 GTGAACGTGGGCTGCGGGCCAGCCGAGGAGCGGGTGTCTGCTGACCGGCCTCCATGCTGTGCGCCGACATC
 CACTGCGAGAACTGCAAGACCACTTTGGGCTGGAATATGAACAGGCCTTTGAGAGCAGCCAGAAGTAC
 AAAGAGGGGAAGTACATCATTGAACTCAACCACATGATCAAAGACAACGGCTGGGACTGA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_031477
Insert Size:	474 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_031477.4</u>
RefSeq Size:	1588 bp
RefSeq ORF:	474 bp
Locus ID:	83719
UniProt ID:	<u>P61236</u>
Cytogenetics:	16p11.2
Domains:	Yippee
MW:	17.5 kDa
Gene Summary:	Involved in proliferation and apoptosis in myeloid precursor cells.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents a longer transcript and longer protein isoform (1).