

Product datasheet for SC317279

REEP5 (NM_005669) Human Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: REEP5 (NM_005669) Human Untagged Clone
Tag: Tag Free
Symbol: REEP5
Synonyms: C5orf18; D5S346; DP1; POB16; TB2; Yip2e; YOP1
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >SC317279 representing NM_005669.
 Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCTGCGGCCATGAGGGAGAGGTTTCGACCGGTTCTGCACGAGAAGAACTGCATGACTGACCTTCTG
GCCAAGCTCGAGGCCAAAACCGGCGTGAACAGGAGCTTCATCGCTCTTGGTGTCATCGACTGGTGGCC
TTGTACCTGGTGTTCGGTTATGGAGCCTCTCTCTGCAACCTGATAGGATTTGGCTACCCAGCCTAC
ATCTCAATTAAGCTATAGAGAGTCCCAACAAAGAAGATGATACCCAGTGGCTGACCTACTGGGTAGTG
TATGGTGTGTTTCAGCATTGCTGAATTCCTCTGATATCTTCTGTGATGTTCCCTTCTACTACATG
CTGAAGTGTGGCTTCTGTTGTGGTGCATGGCCCCGAGCCCTTAAATGGGGCTGAAGTCTCTACAAG
CGCATCATCCGTCTTTCTTCTGAAGCAGAGTCCCAGATGGACAGTGTGGTCAAGGACCTTAAAGAC
AAGGCCAAAGAGACTGCAGATGCCATCACTAAAGAAGCGAAGAAAGCTACCGTGAATTTACTGGGTGAA
GAAAAGAAGAGCACCTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI
Plasmid Map: □
ACCN: NM_005669
Insert Size: 570 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_005669.4</u>
RefSeq Size:	3116 bp
RefSeq ORF:	570 bp
Locus ID:	7905
UniProt ID:	<u>Q00765</u>
Cytogenetics:	5q22.2
Domains:	TB2_DP1_HVA22
Protein Families:	Druggable Genome, Transmembrane
MW:	21.5 kDa
Gene Summary:	May promote functional cell surface expression of olfactory receptors.[UniProtKB/Swiss-Prot Function]