

Product datasheet for SC317307

FKSG24 (MPV17L2) (NM_032683) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FKSG24 (MPV17L2) (NM_032683) Human Untagged Clone
Tag:	Tag Free
Symbol:	FKSG24
Synonyms:	FKSG24
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317307 representing NM_032683. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCGCGGGGTGGCTGGCGCCGGCTACGCCGCTGTTATCCGCGGGGCAGCTTCTATTCCAGGGCCGC
 GCGCTGCTCGTCACTAACACGCTGGGCTGCGCGCGCTCATGGCGGCCGGTGATGGCGTGCGCCAGTCC
 TGGGAGATCCGCGCCCGGCCGCGCCAGGTTTTCGACCCACGGCGCTCCGCGAGCATGTTTGCGGTGGGC
 TGCAGCATGGGTCCCTTCTGCACTACTGGTACTTGTGCTGGACCGCTATTCCTGCGTCTGGCCTC
 CGAGGCTTCCAAATGCTCTCAAGAAGGTCCTCGTGGATCAGCTGGTAGCCTCTCCATTGCTGGGCGTC
 TGGTACTTCTTGGGCCTTGGCTGCCTGGAGGGTCAGACAGTGGGTGAGAGCTGCCAGGAGCTGCGGGAG
 AAGTTCTGGGAATTCTACAAGGCAGACTGGTGCCTGTGGCCTGCTGCGCAGTTCGTGAACCTCCTCTTC
 GTGCCCCCAATTTTCGAGTCACCTACATCAACGGCCTGACGCTGGGCTGGGACAGCTACCTGTCCTAC
 TTGAAGTACCGGAGCCAGTTTCTCTGACACCCCAAGGCTGTGTGGCCCTGGACACCCGAGCAGACTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_032683
Insert Size:	621 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_032683.2
RefSeq Size:	1360 bp
RefSeq ORF:	621 bp
Locus ID:	84769
UniProt ID:	Q567V2
Cytogenetics:	19p13.11
MW:	23.2 kDa
Gene Summary:	Required for the assembly and stability of the mitochondrial ribosome (PubMed:24948607). Is a positive regulator of mitochondrial protein synthesis (PubMed:24948607).[UniProtKB/Swiss-Prot Function]