

Product datasheet for **SC317405**

RBM11 (NM_144770) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RBM11 (NM_144770) Human Untagged Clone
Tag:	Tag Free
Symbol:	RBM11
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_144770, the custom clone sequence may differ by one or more nucleotides ATGTTCCCTGCTCAGGAGGAGGCCGACAGGACCGTGTGGGAAATTTAGAGGCCCGA GTTCCGGGAAGAGATTCTGTACGAGCTGTTCCCTCAGGCGGGGCCACTAACCAAAGTGACT ATATGCAAAGACAGAGAAGGAAAGCCAAAGTCTTTGGATTTGTCTGCTTTAAACACCCA GAATCGGTGTCTTATGCCATAGCTTTGCTGAATGGAATTCGTTTATATGGAAGACCAATT AACGTGCAGTATCGATTTGGGAGTTCTCGCTCTTCTGAACCAGCTAACCAAAGTTTTGAG AGCTGTGTTAAGATAAAATCACACAACACTACAGGAATGAAGAAATGTTGGTGGGCAGATCT TCCTTTCCCATGCAGTATTTTCCAATTAATAACTTCTTTACCTCAAGAATATTTTCTC TTCAGAAAGATGCAGTGGCATGTGTATAATCCAGTGCTGCAGCTTCCTTACTATGAAATG ACAGCTCCACTTCTAATAGTGATCCGTGTCTTCTCACTGAATCATGTTCCAGATCTT GAGGCTGGACCCAGCTCATATAATGGACTCACCAACAACCAAGTGACTCTGACCTTTAT CAGATGACAGCTCCACTTCTAATAGTGATCCGTGTCTTCTCACTGAATCATGTTCCA GATCTTGAGGCTGGACCCAGCTCATATAATGGACTCACCAACAACCAAGTGACTCTGAC CTTTATCAGATGAATAAACGAAAGAGACAAAAGCAAACAAGTGATAGTAGTAGCACA GACAACAACAGAGGCAACGAATGTAGCCAAAAGTCCGAAAGTCTAAGAAGAAGAAAAGA TAC
Restriction Sites:	Please inquire
ACCN:	NM_144770
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_144770.3](#), [NP_658983.3](#)

RefSeq Size: 1997 bp

RefSeq ORF: 846 bp

Locus ID: 54033

UniProt ID: [P57052](#)

Cytogenetics: 21q11.2

Gene Summary: Tissue-specific splicing factor with potential implication in the regulation of alternative splicing during neuron and germ cell differentiation. Antagonizes SRSF1-mediated BCL-X splicing. May affect the choice of alternative 5' splice sites by binding to specific sequences in exons and antagonizing the SR protein SRSF1.[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) encodes isoform 1.