

Product datasheet for SC317479

RFPL2 (NM_006605) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	RFPL2 (NM_006605) Human Untagged Clone
Tag:	Tag Free
Symbol:	RFPL2
Synonyms:	RNF79
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317479 representing NM_006605. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGAAAAGGTTGTCACTTGTCACTAAACAGGCTTTCACCTCATGGAAATTTTTTACCTTGTGTACT
 TTTCCCTGGCAGTGGACATGGCTGCACTCTTCCAAGAAGCAAGCAGCTGTCCGTCTGCTCAGACTAT
 CTGGAAAAACCAATGTCCTGGAGTGTGGATGCGCCGTCTGCCTCAAGTGCATTAATTCAGTGCAGAA
 GAGCCCCATGGGAGGATCTACTTTGCTGTTGCTCTTCCATGGTCTCTCGGAAGAACAAATCAGGCGC
 AATCGGCAGCTAGAGAGGCTGGCTTCCACATCAAGGAAGTGGAGCCCAAGCTGAAGAAGATTCTGCAG
 ATGAACCAAGGATGCGGAAGTTCCAAGTGGATATGACCTTGGATGCCAACACAGCCAACAACCTTCCTC
 CTATTTCTGACGACCTCAGGAGCGTCCGAAGTGGGCGCATCAGACAGAATCGGCAAGACCTTGCCGAG
 AGATTTGACGTGTCCGTTTGCATCCTGGGCTCCCTCGCTTTACCTGTGGCCGCCACTGTGGGAGGTG
 GACGTGGGAACAAGCACAGAATGGGACCTGGGAGTCTGCAGAGAATCTGTTACCCGCAAAGGGAGGATC
 CAGCTGACCACAGAGCTTGGATTCTGGACTGTGAGTTTGAAGGATGGAGGCCGCTCTCTGCCACCACG
 GTGCCGCTGACTTTCCTTTCGTAGACCGCAAGTTACAGCGAGTGGGATTTTTCTGGATATGGGCATG
 CAGAACGTTTCCTTTTGTGCTGAAAGTGGTCCCATGTCTATACATTAGGAGCGTATCTGCTGAG
 GAGCCATTGCGCCATTTTGGCTCCTTCAGTTCCACCTAATGGTGATCAAGGTGTCTTGAGCATCTGT
 CCTTTGATGAAGTCAAGCACTACTGATGCTCCAGTCCGTCCTGGGAGGCCAAATAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_006605
Insert Size:	954 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_006605.3</u>
RefSeq Size:	1496 bp
RefSeq ORF:	954 bp
Locus ID:	10739
UniProt ID:	<u>O75678</u>
Cytogenetics:	22q12.3
Protein Families:	Druggable Genome
MW:	35.5 kDa