

Product datasheet for SC317512

FAM131A (NM_144635) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM131A (NM_144635) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM131A
Synonyms:	C3orf40; FLAT715; PRO1378
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317512 representing NM_144635. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTT TAGTGAACCGTCAGAATTTTGT AATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC GCGATCGCC
ATGTTTCTGGCCACTTTGTCTTTCTCCTCCCCTTGCACATCCCTTTGGAAGTGTTCCTGTGAGTAC
ATGCTGGGGTCTCCCTTTCTCCCTTGCTCAGGTGAATCTCAGCCCCTTCTCCACCCAAAGGTTAC
ATGGATCTAACTACTGCCACCTTCCACCTCCCTGCACCTGTGCTCCCTGGCCTGGTCTTTACCAGG
CTTCTCCACCTCCCTATCTCCAGGTATTTCCAGGTGGTGAAGGACCACGTGACCAAGCCTACCGCC
ATGGCCAGGGCCGAGTGGCTCACCTCATTGAGTGAAGGGCTGGAGCAAGCCGAGTGACTCACCTGCT
GCCCTGGAATCAGCCTTTTCTCCTATTAGACCTCAGCGAGGGCGAACAAGAGGCTCGCTTGCAGCA
GGAGTGGCTGAGCAGTTTCCATCGCGGAAGCCAAGCTCCGAGCATGGTCTTCGGTGGATGGCGAGGAC
TCCACTGATGACTCCTATGATGAGGACTTTGCTGGGGGAATGGACACAGACATGGCTGGGACGTGCC
CTGGGGCCGCACCTCCAGGACCTGTTACCCGGCCACCGGTTCTCCCGCCTGTGCGCCAGGGCTCCGTG
GAGCCTGAGAGCGACTGCTCACAGACCTGTCCCAGACACCTGTGCTTAGTCTGTGCAGCCTGGAG
GATGGGTTGTTGGGCTCCCCGGCCCGGCTGGCCTCCAGCTGCTGGGCGATGAGCTGCTTCTCGCCAAA
CTGCCCCCAGCCGGAAAGTGCCTTCCGACGCTGGGCCACTGGAGGCCAGGACTCACTCTACAAC
CTCTGCCACCACCTAACGGGACGCTGGGAACGGCAGCGGCAAGCCTCTGACCTGGCCTCTTCTGGGGTG
GTGTCCTTAGATGAGGATGAGGCAGAGCCAGAGGAACAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites:	Sgfl-MluI
ACCN:	NM_144635
Insert Size:	1008 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_144635.3
RefSeq Size:	2566 bp
RefSeq ORF:	1008 bp
Locus ID:	131408
UniProt ID:	Q6UXB0
Cytogenetics:	3q27.1
Protein Families:	Secreted Protein
MW:	36.1 kDa