

Product datasheet for **SC317523**

REM2 (NM_173527) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	REM2 (NM_173527) Human Untagged Clone
Tag:	Tag Free
Symbol:	REM2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317523 representing NM_173527. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGCATCGCC**
ATGCACACGGACCTGGACACAGACATGGACATGGACACAGAAACCACAGCACTCTGCCCTCTGGCAGC
CGCCGGGCTCCCTCCAGGGACGCCCACACCAGAAGCAGATGCCACGCTACTAAAGAAGTCAGAGAAA
CTGTTGGCAGAGTTGGACCGAGCGGTTACCTCTGCCCTGGGGCCCCAGACGAAGAGGCAGTATG
CCTGTCCCCTACAAGCACCAGCTCCGGCGGGCCAGGCTGTAGATGAACCTTGACTGGCCACCTCAGGCC
TCATCCTCTGGCTCGTCTGACTCCTTGGGCTCAGGGGAGGCAGCCCTGCTCAAAGGATGGCATCTTC
AAGGTCATGCTAGTGGGGGAGAGCGGCTGGGCAAGAGCACCTAGCAGGCACTTTTGGTGGTCTCCAG
GGAGACAGTGCTCACGAACCGGAGAACCAGAGGATACCTATGAGAGACGCATCATGGTGGATAAAGGAG
GAAGTGACTCTAGTCGTTTATGACATCTGGGAACAGGGGATGCAGGAGGTGGCTGCGGGACCACTGC
CTTCAGACCGGGACGCCTTTCTCATCGTCTTCTCAGTCACCGACCGACGGAGTTTCTCAAAGTTCCA
GAGACCTACTTCGGCTCCGGGCTGGGAGGCCGACACGACCTACCGTTATCCTCGTTGGAAACAAG
AGCGACTTGCCCGCTCCCGGGAGGTATCACTGGAGGAGGGCCGCCACCTGGCCGGGACGCTGAGCTGC
AAGCACATCGAGACGTCGGCCGCACTGCACCACAACACGAGGGAGCTCTCGAGGGCGCGGTGCGCCAG
ATCCGGCTGCGCGGGGGCCGAAACCACGCCGAGGCCAGAGGCCCGATCCGGGCAGCCCCGAGGGCCCT
GCGCCACCTGCACGCCGCGAGAGCCTACCAAGAAAGCCAAGAGGTTCTCGCCAACCTGGTGCCGCGC
AACGCCAAGTTCTTCAAGCAGCGCTCCAGGTCTGTCTACGACCTCTCGGTGCTCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_173527
Insert Size:	1023 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_173527.2</u>
RefSeq Size:	1885 bp
RefSeq ORF:	1023 bp
Locus ID:	161253
UniProt ID:	<u>Q8IYK8</u>
Cytogenetics:	14q11.2
MW:	37.1 kDa
Gene Summary:	Binds GTP saturably and exhibits a low intrinsic rate of GTP hydrolysis.[UniProtKB/Swiss-Prot Function]