

Product datasheet for **SC317534**

ZNF187 (ZSCAN26) (NM_152736) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF187 (ZSCAN26) (NM_152736) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZSCAN26
Synonyms:	SRE-ZBP; SREZBP; ZNF187
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317534 representing NM_152736. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCAACAGCATTGGACCCAGACCAGCCAAAGAAACAAAAAATACTTGTGGAGGAGATGGCCCTCTG
AAAGGAGTACAGGAACAGCAGGTTCCGGCATGAGTGTGAAGTTACAAAGCCTGAGAAAGAGAAGGGTGA
GAGACAAGGATTGAGAATGGGAAGCTTATTGTAGTAACAGACTCTTGTGGAAGAGTAGAGTCATCTGGG
AAAATATCTGAACCCATGGAGGCTCATAATGAGGGCTCTAACTTGAAAGGCATCAGGCCAAGCCAAA
GAGAAGATTGAGTATAAATGCTCAGAACGTGAGCAGAGATTCATCCAGCACTTGGACCTGATTGAACAT
GCGAGTACACACACGGGAAAGAACTCTGCGAGTCTGATGTGTGTCAGAGTTCAGTCTTACAGGACAT
AAGAAAGTCCCTCTAGAGAGAAAGTTCATCAGTGTGATGAGTGTGGAAAGCCTTTCAGAGGAGTTCA
CACCTCGTCAGACATCAGAAAATCCATCTTGGTGAGAAGCCTTATCAGTGCAATGAGTGTGGCAAAGTC
TTTAGCCAGAATGCAGGCCTTTTGGAACTCTCAGAATTCATACTGGAGAGAAACCTTATCTATGTATC
CATTGTGGAAAAATTTAGGCGCAGCTCTCACCTTAATCGACATCAGAGAATTCACAGTCAGGAGGAG
CCCTGTGAGTGAAGGAGTGTGGAAAAACCTTTAGTCAGGCCTTACTCCTCACCCACCATCAGAGAATC
CATAGTCACTCCAAAAGCCATCAATGTAACGAGTGTGGAAAAGCTTTCAGTTTGACCTCAGACCTTATT
CGACACCACAGAATTCATACTGGAGAAAAACCTTTCAAGTGTAAACATATGCCAGAAAGCCTTCCGACTA
AACTCACACCTTGCTCAGCATGTAAGAATCCACAATGAAGAAAAACCTATCAGTGTAGTGAATGTGGA
GAAGCCTTCAGGCAAAGGTCAGGTCTTTTTCAACATCAGAGATATCACCACAAAGACAAACTGGCTTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_152736

Insert Size: 1035 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_152736.5](#)

RefSeq Size: 2445 bp

RefSeq ORF: 1035 bp

Locus ID: 7741

UniProt ID:	<u>Q16670</u>
Cytogenetics:	6p22.1
Domains:	zf-C2H2
Protein Families:	Transcription Factors
MW:	39.7 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (c) uses an alternate in-frame splice site in the 5' coding region, compared to variant a. Variants c and d encode the same isoform (c), which is shorter than isoform a.