

Product datasheet for SC317549

FAM98C (NM_174905) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM98C (NM_174905) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM98C
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL5</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>OriGene sequence for NM_174905 edited

```

GCGCCGAGACCACACTTTCATGGAGGCGGTGAAGCGGAAGCGTGGGAGGGGGCCGCGGT
GGCCCAGGACCTGCTGGCTCTGGGTATGGAGGTGTCCCGGGGGCGGCGTCGCGGGGCGC
CTCATGCCAGACTTCAGGGGGCTGTGCGTGCGGCTGGCGGCGGAGCTGGCGACGCTGGG
CGCCCTCGAGCAGCAGCAGAGGCGGGCGCGGAGGTGCTGAGCGCCGGCGACGGCCCTGG
CGCGGAGGAGGACTTTCGCGGCAGCTCGGCAGCCTGCTGCGGGAGCTGCACTGCCCGGA
TCGCGCGCTCTGCGGCGGGGATGGCGCGGCTGCGCTTCGGAACCCGGTGCCGACTGCG
CCTGCTGCGCTTCTCTGCTCAGAGCTCCAAGCCACCCGCTCCTGTGCCTCCGCTCTCT
GCTGGATCCGAGTCTAGGCCACCCCTTGGTGAAGGGGTAGTGGAGGGAGCCGGCATGGT
CCAAGAAGTGGACCTTACCCTCCAAGCCCTGGGGCTGCCAGACCTGCACAGGGACCCC
CGCCAGCCAGCTGCTGCAGGAGTTGCATGCTAAGATCTCAGAGCTGCAGCCTTCTCTGCC
CCCAGGGTCCCTGCAGCCCTCCTCAGCTGCTCGCTAGATGCACCCAGATGGGAAGCGTT
GGAGTCTCTGTCCAAAGCCTCAGAGATCAGTACCGCTGCCGCCGCTGCCTCCTCCTCAA
GCGCCTTGACCTCACTACATCTGCTTCCACTGGAGTGACCGGGCAGAGGCCAAGGAGA
GGCCATGAGGGCAGTGCTGATCCCAATTCGAGAGGTTCTGACCCAGAATCGGACATCTC
CATTGCACACGTTCTGGCTGCCGAGCCGACCTGTCTGTCTCGTCCCAGCCACCAGCGT
GGCTGTCCGAGAGGGACCTGCTGTGCCATCAACAAGGTGCTTATGGGCAACGTTCCAGA
CCGGGGGGGGCCGCCAAATGAGCTGGAGCCTCCATGCCACCTGGAGGAGCCGAAGAGA
GGATGGAGGCCCCAGTGTGGGGTCGCAAGAAGAAGAAGAAGTAAAGGGGGACTGG
TGGTCCGGGGCGGAGGTCCTCATGAGATGCTAATGCTATCGGTAACTCTCTGGGAGTC
CGTTTAAGGCTTTCCTTGGTATTTGGCACCCAAGCCCCCATCTCCTATTCTGAGTCCC
CAAATGCCCTGCTCCATTACGTCAA
  
```

Restriction Sites:	Please inquire
ACCN:	NM_174905
Insert Size:	1200 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	The ORF of this clone has been fully sequenced and found to be a perfect match to NM_174905.3.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_174905.3</u> , <u>NP_777565.3</u>
RefSeq Size:	1291 bp
RefSeq ORF:	1050 bp
Locus ID:	147965
UniProt ID:	<u>Q17RN3</u>
Cytogenetics:	19q13.2