

Product datasheet for SC317561

SH3GL3 (NM_003027) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	SH3GL3 (NM_003027) Human Untagged Clone
Tag:	Tag Free
Symbol:	SH3GL3
Synonyms:	CNSA3; EEN-B2; HsT19371; SH3D2C; SH3P13
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317561 representing NM_003027. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTCGGTGGCCGGGCTGAAAGAAGCAGTTCACAAAGCCAGCCAGCTATTTAGTGAAAAATAAGTGGT
GCTGAAGGAACTAACTAGACGATGAATTTCTTGACATGGAAAGGAAAATAGATGTTACCAATAAAGTT
GTTGCAGAAATCTTTCAAAAACCACTGAATATCTTCAGCCAAATCCAGCATAACAGAGCTAAGCTAGGA
ATGCTGAACACTGTGTCGAAGATCCGAGGGCAGGTGAAGACCACAGGATACCCGACAGCGGAAGGCTTG
CTGGGGGACTGTATGCTGAAATACGGGAAGGAGCTCGGGGAAGACTCCACCTTTGGCAATGCATTGATA
GAAGTTGGTGAATCCATGAAGCTAATGGCTGAGGTGAAAGACTCTTTGATATTAATGTAAGCAAACCT
TTTATTGATCCACTTCAGTTACTACAAGATAAAGATTTAAAAGAGATCGGGCATCACCTGAAAAAGCTG
GAAGGCCGCCGCTGGATTACGATTATAAAAAGAAACGAGTAGGTAAGATACCAGACGAAGAAGTCAGA
CAAGCGGTAGAAAAATTTGAAGAGTCAAAGGAGTTGGCTGAAAGAAGCATGTTTAACTTTTTAGAAAAAT
GATGTAGAACAAGTCAGCCAGTTGGCTGTGTTTCATAGAGGCAGCATTAGACTATCACAGACAGTCCACA
GAGATTCTGCAGGAGCTGCAGAGCAAGCTACAGATGCGAATATCAGCTGCATCCAGTGTCCCCAGACGA
GAATAACAAGCCAAGGCCGTGAAAAGGAGTTCTAGTGAGCTCAATGGAGTTTCCACCACCTCTGTAGTG
AAGACGACAGGTTCTAACATTCATGGACCAGCCCTGCTGTGCGTGTCTCTATGACTTTGAGCCAGAA
AACCAAGGAGAAATTAGGATTTAAAGAAGGGGACATCATTACATTAACCAATCAAATAGATGAAAACCTGG
TATGAAGGAATGATACACGGAGAATCGGGATTCCTCCCATTAATTTACGTGGAAGTGATCGTGCCTTTA
CCTCAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_003027

Insert Size: 1044 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_003027.4](#)

RefSeq Size: 1843 bp

RefSeq ORF: 1044 bp

Locus ID: 6457

UniProt ID:	<u>Q99963</u>
Cytogenetics:	15q25.2
Domains:	SH3, BAR
Protein Pathways:	Endocytosis
MW:	39.3 kDa
Gene Summary:	Implicated in endocytosis. May recruit other proteins to membranes with high curvature (By similarity).[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (EEN-B2-L1) encodes isoform a.