

Product datasheet for **SC317575**

FAM102B (NM_001010883) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	FAM102B (NM_001010883) Human Untagged Clone
Tag:	Tag Free
Symbol:	FAM102B
Synonyms:	SYM-3B
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317575 representing NM_001010883. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
ATGATGAAGAAGAAGAAGTTTAAGTTTAAGGTGGACTTCGAGCTCGAGGAGCTCTCCTCAGTGCCCTTC
GTCAATGGGGTCTCTTTGCAAGATGCGGCTGCTGGACGGCGGCAGCTTCACCGCCGAGTCCTCCAGG
GAGGTGGTACAAGCAAAGTGTGTCGCTGGAGAAAGAAGTTCTCATTTATGTGCAAAATGAGTGCAAGT
GCTGCCACAGGCATCCTAGATCCTTGTATCTACAGAGTATCCGTGAGGAAGGAATTAAGGTGGAAAA
GCTTATGCAAGCTGGGCTTTGCAGATCTAAACCTGGCAGAGTTTGTGGATCAGGAAATACCACTCGC
CGCTGTTTACTGGAAGGCTATGATACCAAAAATACAAGACAGGATAATTCCATTCTTAAAGTTTGTATC
AGTATGCAACTGATGTCTGGTGACCCATGTTTTAAACGCCTCCCTCCACTTCAATGTCTATACCAATT
GCTGGTGAATCTGAATCTTTGCAAGAAGATAGAAAAGGTGGAGAGACTCTCAAAGTGCATCTTGAATA
GCAGATCTTTAGCAAGAGTGCCTCTGTTCCAGACGAACTTGGTGCCTGTGGACATTCTAGAATCA
AGCTATGCAAGTCAGCAGTCAAAAGTATCAGGGTATAGCACGTGTCACTCCCGGTCATCTAGTTTCTCT
GAGCTCTGCCACAGGAGAAATACCTCAGTGGGAAGCACATCAACAGGAGTTGAAAGTATTCTAGAGCCA
TGATGATGAAATTGAGCAGAAAATCGCTGAGCCAAATCTTGATACAGCTGATAAAGAAGATACAGCTTCA
GAAAACTCAGCAGATGCCCAGTGAAACAAGATTCTGTAGAATCTCAGCTGAAGCGAGTTGATGACACC
AGGGTGGATGCAGATGACATTGTAGAGAAAATATTACAAAGTCAAGACTTCAGCCTAGATTCCAGTGCG
GAAGAAGAAGGACTAAGGCTATTTGTGGTCTCGGGGAAGTACAACCTTTGGCAGTCATCATCTTCCA
AATAGGTGGGTCTGGAGCTTATGAACAAGTTGTGATCAAACGC**TAG**
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC

Restriction Sites:	SgfI-MluI
ACCN:	NM_001010883



[View online »](#)

Insert Size:	1083 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001010883.2</u>
RefSeq Size:	5359 bp
RefSeq ORF:	1083 bp
Locus ID:	284611
UniProt ID:	<u>Q5T8I3</u>
Cytogenetics:	1p13.3
MW:	39.3 kDa