

Product datasheet for **SC317578**

ATAD1 (NM_032810) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ATAD1 (NM_032810) Human Untagged Clone
Tag:	Tag Free
Symbol:	ATAD1
Synonyms:	AFDC1; FNP001; HKPX4; Msp1; THORASE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317578 representing NM_032810. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTACATGCTGAAGCCTTTTCTCGTCTTTGAGTCGGAATGAAGTTGTTGGTTAATTTTCCGTTTG
ACAATATTTGGTGCAGTGACATACTTTACTATCAAATGGATGGTAGATGCAATTGATCCAACCAGAAAG
CAAAAAGTAGAAGCTCAGAAACAGGCAGAAAACTAATGAAGCAAATGGAGTAAAAATGTGAAGCTC
TCAGAATATGAAATGAGTATTGCTGCTCATCTTGTAGACCCTCTAATATGCATGTTACTTGGAGTGAT
ATAGCAGGTTTAGATGATGTCATTACGGATCTGAAAGACACAGTCATCTTACCTATCAAAAAGAAACAT
TTGTTTGAGAATCCAGGCTTCTGCAGCCTCCAAAAGGTGTTCTTCTCTATGGGCCTCCAGGCTGTGGT
AAAACGTTGATTGCCAAGGCCACAGCCAAAGAAGCAGGCTGTCGATTTATTAACCTTCAGCCTTCGACA
CTGACCGATAAGTGGTATGGAGAATCTCAGAAATTTGGCTGCTGCTGTCTTCTCCCTTGCCATAAAGCTA
CAACCATCCATCATCTTTATAGATGAAATAGACTCCTTTCTACGAAACCGTTCAAGTTCTGACCATGAA
GCTACAGCCATGATGAAAGCTCAGTTTATGAGTCTCTGGGATGGATTGGATACTGATCACAGCTGCCAG
GTCATAGTAATGGGAGCTACCAATCGTCCTCAGGACCTTGACTCGGCTATAATGAGAAGAATGCCTACA
AGATTTTCATATCAACCAGCCTGCTTTAAAACAGAGAGAAGCAATCCTGAAACTCATCTTGAAAAATGAA
AATGTGGATAGGCATGTAGACCTGCTAGAAGTTGCCAGGAACTGATGGGTTTTTCAGGAAGTGACCTA
AAAGAGATGTGTCGAGATGCTGCCCTCCTGTGTTAGAGAATATGTTAATTCTACATCAGAAGAAAGC
CATGACGAAGATGAAATTCGGCCTGTTCAACAGCAGGACCTGCATCGGGCAATTGAAAAGATGAAGAAA
TCAAAGGATGCAGCATTTTCAAGATGTTTTAACACATGTTTGTGTTAGATTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_032810

Insert Size: 1086 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_032810.3](#)

RefSeq Size: 4711 bp

RefSeq ORF: 1086 bp

Locus ID: 84896

UniProt ID: [Q8NBU5](#)

Cytogenetics: 10q23.31

MW: 40.7 kDa

Gene Summary: ATPase that plays a critical role in regulating the surface expression of AMPA receptors (AMPA), thereby regulating synaptic plasticity and learning and memory. Required for NMDA-stimulated AMPAR internalization and inhibition of GRIA1 and GRIA2 recycling back to the plasma membrane; these activities are ATPase-dependent (By similarity).
[UniProtKB/Swiss-Prot Function]