

## Product datasheet for SC317581

### NDRG3 (NM\_022477) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	NDRG3 (NM_022477) Human Untagged Clone
Tag:	Tag Free
Symbol:	NDRG3
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317581 representing NM_022477. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGATGAACTTCAGGATGTTTCAGCTCACAGAGATCAAACCCTTCTAAATGATAAGGAACATGATATA
GAAACAACCTCATGGTGTGGTCCACGTCCTATAAGAGGCTTACCCAAAGGAAACAGACCAGTTATACTA
ACATATCATGACATTGGCCTCAACCATAAATCCTGTTTCAATGCATTCTTTAACTTTGAGGATATGCAA
GAGATCACCCAGCACTTTGCTGTCTGTGCATGTGGATGCCCCAGGCCAGCAGGAAGGTGCACCCTCTTTC
CCAACAGGGTATCAGTACCCACAATGGATGAGCTGGCTGAAATGCTGCCTCCTGTTCTTACCCACCTA
AGCCTGAAAAGCATCATTGGAATTGGAGTTGGAGCTGGAGCTTACATCCTCAGCAGATTTGCACTCAAC
CATCCAGAGCTTGTGGAAGGCCTTGTGCTCATTAAATGTTGACCCTTGCCTAAAGGCTGGATTGACTGG
GCAGTTCCAAACTCTCTGGCCTGACAACCAATGTTGTGGACATTATTTGGCTCATCACTTTGGGCAG
GAAGAGTTACAGGCCAACCTGGACCTGATCCAAACCTACAGAATGCATATTGCCAAAGACATCAACCAA
GACAACCTGCAGCTCTTCTTGAATTCCTACAATGGACGCAGAGACCTGGAGATCGAAAGACCCATACTG
GGCCAAAATGATAACAAATCAAAAACATTAAGTGTCTACTTTACTGGTGGTAGGGGACAATTCGCCT
GCAGTTGAGGCTGTGGTGAATGCAATCCCGCCTGAACCCTATAAATACAACCTTGCTAAAGATGGCG
GACTGTGGGGGACTGCCCGAGTGTTCAGCCTGGGAAGCTCACCGAGGCCTTCAAGTACTTTTTGCAG
GGAATGGGCTACATACCATCTGCCAGCATGACTCGGCTCGCCGATCACGAACCCACTCAACCTCGAGT
AGCCTCGGCTCTGGAGAAAGTCCCTTACGCCGGTCTGTACCAGCAATCAGTCAGATGGAACCAAGAA
TCCTGTGAGTCCCTGATGTCCTGGACAGACACAGACCATGGAGGTGTCCTGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_022477

**Insert Size:** 1092 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_022477.3](#)

**RefSeq Size:** 2995 bp

**RefSeq ORF:** 1092 bp

Locus ID: 57446  
UniProt ID: [Q9UGV2](#)  
Cytogenetics: 20q11.23  
Domains: Ndr  
MW: 40 kDa