

Product datasheet for **SC317595**

ZNF707 (NM_173831) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	ZNF707 (NM_173831) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF707
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317595 representing NM_173831. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGACATGGCCAGGAGCCAGTGACCTTCAGGGACGTGGCCATCTACTTCTCAAGGGAGGAGTGGGCG
TGTCTGGAACCCAGCCAGAGGGCCCTCTACCGGACGTGATGCTGGACAACCTCAGCAGTGTGGCTGCT
CTGGGATTTTGCAGCCCCAGACCAGACCTCGTCTCTCGCTGGAACAGTGGGAGGAGCCGTGGGTTGAA
GACCGGGAGAGACCTGAGTTCAGGCAGTGCAGAGGGGACCCCGCCAGGGGCAAGGAAGTCTGCAGAC
CCCAAGAGACCTTGTGATCATCCAGCTTGGGCTCACAAGAAAACCCACGTGCGGCGAGAAAGAGCCAGG
GAAGGAAGCAGCTTAGGAAGGGCTTCAGGCTGGACCGGATGACGGGCAGCTTCCAGAGCTGCTCCA
GAAAGGACAGACGCCAAGCCACGGCTTCCCGTGTGAGGTGCTCACGCAGCGTTGTGGCGGCGGCCG
GGCCGAGAGAGCGCCGGAAGCAGCGCGCAGTAGAGCTGTCATTCTGCGGCACGTGCGGGAAGGCG
CTCAGCTGCCACAGCCGCTGCTCGCTCACCAGACGGTGCACACGGGAACCAAGGCTTCGAGTGGCCC
GAGTGGCGCCAGACCTTCCGGTGGGCTTCAAACCTGCAGCGCCACCAGAAGAACCACGCGCGAGAAG
CCCTTCTGCTGCGAGGCCTGCGGGCAGGCGTTCAGCCTGAAGGACCGCTGGCTCAGCACCAGGTC
CACACCGAGCACAGGCCCTACTCGTGTGGCAGTGTGGGAAAGCCTTCAAGCAGAAGTCCAACCTTCTC
AGACACCAGCTGGTGCACACCGGGAGCGGCCGTTCTACTGCGGACTGCGGCAAGCCTTCCGGACC
AAGGAGAACCTCAGCCACCAGAGGGTCCACAGCGGGGAGAAGCCCTACACCTGTGCCGAGTGGCGC
AAGTCCTTCCGGTGGCCCAAGGGCTTCAGCATCCACCGGAGGCTGCACCTGACGAAGAGGTTCTACGAG
TGGCGCCACTGTGGGAAAGGCTTCCGTCACCTGGGTTCTTACGCGGCATCAGAGGACTCACAGGCAC
GGGAGGTGTAG
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGCGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_173831

Insert Size: 1116 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173831.3](#)

RefSeq Size: 2302 bp

RefSeq ORF: 1116 bp

Locus ID: 286075

UniProt ID:	<u>Q96C28</u>
Cytogenetics:	8q24.3
Protein Families:	Transcription Factors
MW:	43.1 kDa
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) differs in the 5' UTR, compared to variant 4. Variants 1, 2, 3, 4, and 5 encode the same protein (isoform a).