

## Product datasheet for SC317601

### WDR86 (NM\_198285) Human Untagged Clone

#### Product data:


Product Type:	Expression Plasmids
Product Name:	WDR86 (NM_198285) Human Untagged Clone
Tag:	Tag Free
Symbol:	WDR86
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317601 representing NM_198285. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGGGGGCGGGGGTCCGGCCCTGAGGGTCTGCGCCGACCAACCGGGGGGCATCAACTGGCTGAGCCTG
AGCCCCGACGGGCAGCGCCTGCTGACGGGCAGCGAGGACGGCACGGCCCGGCTCTGGAGCACCGCGGAC
GGCCAGTGCTGCGCGCTCCTGCAAGGACATGAAAGCTATGTGACCTTCTGCCAGCTGGAGGATGAGGCT
GCCTTCACATGCAGCGCCGACTGCACCATCAGGAGGTGGGACGTGCTGACCGGGCAGTGTCTGCAGGTG
TACCGAGGACACACGTCATCGTGAACAGGATCCTGGTTGCCAACACCAGCTCTTCAGCAGCTCCTAT
GACCGGACAGCTCGGGTCTGGAGTGTGGACAAGGGGCAGATGTCCCGGGAGTTCGGGGCCACCGCAAC
TGCGTGCTGACCCCTAGCCTACTCTGCCCGTGGGACCTCCCGAGCACTCCCTGCGCGGAGGAGCCGCG
GCCGGGGGCTTCTGGTGACCGGCAGCACAGATGGCACAGCCAAGGTGTGGCAGGTGGCCAGCGGTGC
TGCCACCAGACGCTGCGGGGCCACCGGTGCACTGCTGTGCCTAGTGCTAGACACGCCCCGCCACACG
GCCTTCACAGGCAGCACCGACGCCACCATCCGTGCCTGGGACATCCTGAGTGGGAGCAGCTGCGGGTG
TTCCGGGAGCACCGGGCTCCGTCATCTGTCTGGAGCTGGTGAACCGACTCGTGTACTCTGGCAGCGCG
GACAGGACCGTCAAGTGTGGCTGGCAGACACAGGGGAGTGTGTGCGCACGTTACGCCCCACAGACGC
AACGTGAGCGCCCTCAAGTACCACGCGGGCACCTTGTTACGGGCAGCGGGACGCTTGCGCCGGGCC
TTCGACGCGCAGTCTGGAGAGCTGCGGAGGGTGTCCGGGGCCACACATTCATCATCAACTGCATCCAG
GTGCACGGCCAGGTGCTCTACACCGCCTCGCACGACGGCCCTCGCCCTCTGGGACGTGCGCGGGCTC
CGAGGTGCCCCGCGGCCCTCCGCCATGCGCAGCCTCTCGCGGCTCTCAGCAACAAGGTGGGCTGC
GCCGCCGCGCCCTGCAGCCGGCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

**Plasmid Map:**


**ACCN:** NM\_198285

**Insert Size:** 1131 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_198285.2](#)

**RefSeq Size:** 2056 bp

**RefSeq ORF:** 1131 bp

**Locus ID:** 349136

UniProt ID: Q86TI4  
Cytogenetics: 7q36.1  
MW: 40.7 kDa