

Product datasheet for **SC317605**

METTL2A (NM_181725) Human Untagged Clone

Product data:


Product Type:	Expression Plasmids
Product Name:	METTL2A (NM_181725) Human Untagged Clone
Tag:	Tag Free
Symbol:	METTL2A
Synonyms:	METTL2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317605 representing NM_181725. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```
GCTCGTTTAGTGAACCGTCAGAATTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCGGCTCCTACCCTGAAGGTGCACCTGCAGTCTCGCCGATAAGAGGCAGCAGTTCGGAAGCCGG
TTCCTGAGAGATCCGGCGCGCGTCTCCACCACAATGCCTGGGACAATGTGGAGTGGTCGGAAGAGCAA
GCCCGCGCGGGAGAGAAAAAGTCCAGGAGAACAGTATCCAGCGGGTGTGCCAGGAGAAACAAGTTGAT
TATGAGATCAATGCCACAAATACTGGAATGACTTCTACAAAATCCACGAAAATGGGTTTTCAAGGAT
AGACATTGGCTTTTTACCGAATTCCTGAGCTGGCACCTAGCCAAAATCAAATCATTGAAGGACTGG
TTCTTGGAGAACAAGAGTGAAGTACCTGAATGTAGAAACAATGAGGATGGACCTGGTTTAATAATGGAA
GAACAGCACAAGTGTCTCAAAGAGCCTTGAACATAAAACACAGACACTTCTGTGGAGGAGAATGTA
ACTCAGAAAATTAGTGACCTGGAATTTGTGCTGATGAGTTTCTGGATCCTCAGCCACCTACCGAATA
CTGGAGGTTGGCTGTGGTGTGGGAAACACAGTCTTTCCAATTTTACAAACGAACAATGACCCAGGACTC
TTTGTATTGCTGTGATTTTTCTCCACAGCTATAGAAGTGGTCCAGACAAATTCAGAATATGATCCT
TCTCGGTGTTTTGCCTTTGTTCCAGACCTGTGTGATGAAGAGAAGAGTTACCCAGTGCCCAAGGGCAGT
CTTGATATTATCATTCTCATATTTGTTCTTTCCAGCAATGTTCCAGACAAGATGCAGAAGGCTATCAAC
AGGCTGAGCAGGCTTCTGAAACCTGGCGGGATGATGCTTCTGCGAGATTACGGCCGCTATGACATGGCT
CAGCTTCGGTTTTAAAAAGGTCAAGTGTCTATCTGGAAATTTCTACGTGAGAGGTGATGGAACAGAGTT
TACTTCTCACACAAGAGGAAGTGGACACGCTTTTACCACCTGCTGGACTGGAAAAAGTTTCAAGACCTG
GTGGATCGCCGACTGCAGGTGAACCGAGGAAAGCAACTGACAATGTACCGGGTTTGGATTAGTGCAAA
TACTGCAAGCCCCTTCTGTCCAGCACCAGTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:


ACCN: NM_181725

Insert Size: 1137 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_181725.3](#)

RefSeq Size: 2519 bp

RefSeq ORF: 1137 bp

Locus ID: 339175

UniProt ID: [Q96I76](#)

Cytogenetics: 17q23.2

Protein Families: Druggable Genome

MW: 43.5 kDa

Gene Summary: Probable S-adenosyl-L-methionine-dependent methyltransferase that mediates 3-methylcytidine modification of some tRNAs.[UniProtKB/Swiss-Prot Function]