

Product datasheet for **SC317623**

PTCD2 (NM_024754) Human Untagged Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | PTCD2 (NM_024754) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | PTCD2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC317623 representing NM_024754. Blue=Insert sequence Red=Cloning site Green=Tag(s) |


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGTCCGAGACAGTATGGCTGCTGCATTTTCGGCCCTCGAATCGAGTTCTCCTGCAGGCGCTGCAGATT
TTGGTGTATCCTGGGGTGGGAGGCTCCGGCTCTGTGAGCTGCCGCTGCCCTCTCGGAGCTAAAAGATAC
CTACTTACAGATAATGTGGTGAATTTAAAAGAATTTCAACAAAAGAAAGTGGCTGTTGCATGTAATCTT
TCTGGCACTAAAGAAACGTATTTTAGAACTTGAAAAAGAACTGACCCAGAACAAGCTCATCTTGAAG
GGGGAGTTGATAACCTTACTACATTTGTGTGAGTCTCGGGACCATGTGGAAGTGGCTAAAAATGTCATT
TACAGGTACCATGCAGAGAACAAAAATTTCACTTTGGGGGAGTATAAATTTGGACCGCTTTTGTGAGG
TTGTGTTACGAGTTGGATCTCGAGGAATCTGCAGTGGAGCTCATGAAAGACCAGCATTTACGAGGTTTC
TTCTCAGACTCCACATCATTCAATATTTTATGATGATATGTTATTTATCAAAGGCAAAATATAAAAGTGT
TTGCAAGTATTGATAGAGATGAAAAACCAAGATGTGAAGTTCACCAAAGATACCTATGTTCTTGCTTTT
GCAATTTGCTACAACTGAATAGCCCTGAGTCTTTCAAATCTGTACTACATTAAGAGAAGAAGCTCTA
CTCAAAGGAGAAATCTCTCCAGGAGAGCATCCTGTTTCGCTGTGGCATTAGCTCTGAATCAGAATGAG
ATGGCAAAGCTGTGTCCATTTTTCTCAAATCATGAATCCAGAAAGCATAGCCTGCATTAATTTAAAT
ATTATAATCCATATCCAGTCAAATATGTTGAAAACCTGATAAAGACTCTAAAAATGCTGCAGAAGGA
AATTTATCAAATTTGTGAAAAGACATGTGTTCTCGGAGGAAGTCTGGCCAAAGTGAAGGAAAAAGTG
AAGGATGTGCCTGCCCTTGTGGCCAAATTTGATGAGATCTATGGGACACTGCACATCACTGGCCAGGTC
ACCACTGATTCTTTGGATGCTGTGCTCTGCCACACCCCGAGGACAGGAAATCTCACACGTTGCTATTA
AACAAAGAGGATGGTCAGCCGTCGCACCTTCCAGCCACTCAGCCAGTCCCTGTTGGCTGAGTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_024754

Insert Size: 1167 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024754.4](#)

RefSeq Size: 2128 bp

RefSeq ORF: 1167 bp

Locus ID: 79810

UniProt ID: [Q8WV60](#)

Cytogenetics: 5q13.2

Domains: PPR

MW: 44 kDa

Gene Summary: Involved in mitochondrial RNA maturation and mitochondrial respiratory chain function.
[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).