

Product datasheet for SC317626

FAM80A (RIMKLA) (NM_173642) Human Untagged Clone

Product data:


| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | FAM80A (RIMKLA) (NM_173642) Human Untagged Clone |
| Tag: | Tag Free |
| Symbol: | RIMKLA |
| Synonyms: | FAM80A; NAAGS; NAAGS-II; NAAGS2 |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >SC317626 representing NM_173642. Blue=Insert sequence Red=Cloning site Green=Tag(s) |

```
GCTCGTTTGTGAACCGTCAGAATTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGTGCTCCAGCTCTGGTTCTGACGGACCGGCATCCCGGAGGACTACCCGAGGTGCAGATCCTG
CGCGCCCTCCGGCAGCGCTGCTCCGAGCAGGACGTGCGCTTCCGGGCGGTGCTTATGGACCAGATCGCC
GTCACCATCGTCGGCGGCCACCTCGGCCCCAGCTAAACCAGAAGGCCCTCACCCTTCCCGGATGTG
GTGCTTGTACGGGTACCCACACCCTCAGTGCAGTCAAGACAGTGCATCACTGTCTGCGACACCTGGAG
AAGCTGGGCTGCCGTTGGTCAATCGCCACAGAGCATCTAAATGCATCAACAATTCTGGACGTTTCA
CAAGAAGTGGCTGGACATGGGGTCCCCATGCCAGACACCTTCTCCTATGGTGGGCATGAAGACTTTTCA
AAAATGATTGATGAAGCTGAGCCCTGGGCTACCCAGTCGTGGTGAAGAGCACAGAGGCCACCGGGGA
AAAGCTGTTTTTCTGGCAAGAGATAAACATCACCTCTCTGACATCTGCCATCTGATCCGCCACGATGTG
CCCTACCTGTTCCAGAAGTACGTGAAGGAGTCCCATGGAAGGACATCCGGTGGTGGTGGTAGGGGGC
CAGGTCATAGGCTCTATGCTTCCGCTGCTCCACTGATGGACGGATGCAGAGCAACTGCTCTCTCGGTGGC
GTGGGCGTCAAGTGTCCGCTGACAGAACAAGGCAAGCAGTTGGCTATTCAGGTGTCCAACATCCTAGGC
ATGGACTTCTGTGGCATTGATCTCCTTATCATGGACGATGGCTCCTTTGTGGTGTGTGAGGCAAATGCT
AATGTTGGCTTCTAGCCTTTGACCAGGCATGCAACTTAGATGTGGTGGGATCATTGCAGACTATACC
ATGTCCTTGCTGCCAAATAGGCAGACTGGAAGATGGCTGTCTCCAGGACTGTCGAGTCCAAGGGGAG
AAGAACGAGCCGGATGGCTGTGCTTCAAGCTCAGGGAGTTGCAGAGAGCGTCTATACCATCAACAGTGGG
TCTACCTCTAGTAAAAGTGAAGCTGAACTGGGAGAGATCCGGGATTCCTCAGCAAGCACAATGGGGGCC
CCACCCTCCATGCTGCCGAACCTGGCTACAACATTAACAACAGGATTGCTTCTGAGTTAAAACCTAAG
TGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: Sgfl-MluI



[View online >](#)

Plasmid Map:


ACCN: NM_173642

Insert Size: 1176 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173642.3](#)

RefSeq Size: 10583 bp

RefSeq ORF: 1176 bp

Locus ID: 284716

UniProt ID: [Q8IXN7](#)

Cytogenetics: 1p34.2

MW: 42.9 kDa

Gene Summary: Catalyzes the synthesis of N-acetyl-L-aspartyl-L-glutamate (NAAG) and N-acetyl-L-aspartyl-L-glutamyl-L-glutamate.[UniProtKB/Swiss-Prot Function]