

Product datasheet for SC317634

DNAJA4 (NM_018602) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	DNAJA4 (NM_018602) Human Untagged Clone
Tag:	Tag Free
Symbol:	DNAJA4
Synonyms:	MST104; MSTP104; PRO1472
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317634 representing NM_018602. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCCCGGGGGCGCAGTCAGAGCTGGAGCTCCGGGGAATCAGACGGGCAGCCAAAGGAGCAGACGCC
GAGAAGCCAGACACAAGATGGTGAAGGAGACCCAGTACTATGACATCCTGGCGTGAAAGCCAGCGCG
TCCCCGGAGGAGATCAAGAAGGCCTATCGGAAGCTGGCGCTCAAGTACCACCCGGACAAGAACC GGAT
GAGGGCGAGAAGTTAACTCATATCCAGGCATATGAAGTGCTTCAGATCCAAAGAAAAGGGATGTT
TATGACCAAGGCGGAGAGCAGGCAATTAAGAAGGAGGCTCAGGCAGCCCGAGCTTCTTCCACCCATG
GACATCTTTGACATGTTCTTTGGTGGTGGACGGATGGCTAGAGAGAGAAGAGGCAAGAATGTTGTA
CACCAGTTATCTGTAACCTTGAAGATCTATATAATGGAGTCACGAAGAAATTGGCCCTCCAGAAAAT
GTAATTTGTGAGAAATGTGAAGGTGTTGGTGGGAAGAAGGGATCGGTGGAGAAAGTCCCGCTGTGCAAG
GGGCGGGGGATGCAGATCCACATCCAGCAGATCCGGCCGGGCATGGTACAGCAGATCCAGACCGTGTGC
ATCGAGTGCAAGGGCCAGGGTGAGCGCATCAACCCCAAGGACCGCTGCGAGAGCTGCAGCGGGGCCAAG
GTGATCCGTGAGAAGAAGATTATCGAGGTACATGTTGAAAAGGTATGAAAGATGGGCAAAAGATACTA
TTTCATGGAGAAGGAGATCAGGAGCCTGAGCTGGAGCCTGGTATGTCATAATTGTGCTTGTATCAGAA
GATCATAGTGTCTTTAGAGACGAGGCCATGACTTGATCATGAAAATGAAAATTCAGCTTTCTGAAGCT
CTTTGTGGCTTCAAGAAGACGATAAAAACATTGGACAATCGAATCTTTGTTATTACATCCAAAGCAGGT
GAGGTGATAAAGCACGGGGACCTGAGATGCGTGCGGATGAAGGAATGCCCATCTACAAAGCACCCCTG
GAAAAAGGGATTCTGATCATACAGTTTTAGTAATCTTTCTGAAAAAAGTGGCTTTCTCTGAAAAAG
CTTCTCAGTGGAAGCTTTACTCCCTCTCGACAGAAAGTGAGGATTACAGATGACATGGATCAGGTG
GAGCTGAAGGAGTTTTGTCCCAATGAGCAGAACTGGCGTCAGCACAGGGAGGCTACGAGGAGGACGAA
GACGGGCCCCAGGCTGGAGTGCAAGTCCAGACGGCATGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_018602

Insert Size: 1281 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_018602.3](#)

RefSeq Size: 3203 bp

RefSeq ORF: 1281 bp

Locus ID: 55466

UniProt ID: [Q8WW22](#)
Cytogenetics: 15q25.1
Domains: DnaJ_CXXCXGXG
MW: 48 kDa