

Product datasheet for SC317643

PHF23 (NM_024297) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PHF23 (NM_024297) Human Untagged Clone
Tag:	Tag Free
Symbol:	PHF23
Synonyms:	hJUNE-1b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317643 representing NM_024297. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCTGGAAGCCATGGCGGAGCCAGTCCCGAAGATCCACCTCCGACCCTTAAGCCAGAGACTCAGCCA
CCAGAGAAACGGCGGAGAACAATTGAGGATTTCAACAAATTCGCAGTTTGTGTTTGGCATATGCTGGT
TACATTCCCCTAGCAAAGAGGAAAGTACTGGCCAGCCTCTGGCTCCAGCTCTCCATTGCGAGGAGAG
AGTGGCGCCGACAGTGATGGCTGGGACTCGGCCCTCAGATCTTGAACCATCCAGACTTTTGTAAAG
AAAGCAAAGTCATCCAAGAGAAGGGCAGCTCAAGCAGTCCCACCCAGCCAGGACCCCAAGTCCACT
TTCTCTGCTCTGCAGGCCCCGACAGTGCTACCTTGCTTGAGAAGATGAAGCTCAAGGACTCTCTTTT
GATCTGGATGGGCCAAAGTGGCATCTCCTTTGTCCCCACATCCCTGACACATACCTCCCGGCCCCCT
GCTGCTTTACCCCGTGCCTTTCCAGGGGACCTCTCCATCCTCCTCGAAAGAAGGACCGAAAG
AACCGAAAGTTGGGGCCAGGAGCTGGGGCTGGCTTTGGGGTGTTCGGAGGCCTCGGCCAACTCCTGGG
GATGGGAAAAGAGATCTCGAATCAAGAAGAGCAAGAAGCGGAAGTTAAAAAGGCAGAACGGGGGAT
AGACTCCCACCTCCTGGGCTCCCGAGGACCCCCAGTGATACAGACTCTGAAGAGGAGGAGGAAGAG
GAGGAAGAGGAAGAAGAAGAAGATGGCAACAGTGGTAGGGGTGAAGCCCAAGTCCCTGTGCTGCCA
ACACCCCTGAGGCTCCTAGGCCCTGCCACAGTGCACCCTGAAGGAGTCCCTCCTGCTGACAGTGAA
AGCAAGGAGGTGGGCAGCACTGAAACAAGCCAAGATGGAGATGCCAGCTCCAGTGAAGGCGAGATGCGG
GTCATGGACGAGGACATCATGGTAGAATCAGGTGATGACTCATGGGATCTGATCACATGTTACTGTGCGA
AAGCCCTTTGCAGGGCGGCCATGATTGAGTGCAGCCTGTGTGGGACGTGGATCCACCTCTCCTGTGCT
AAGATTAAGAAGACCAACGTCCCGACTCTTTTATTGCCAGAAATGCAAGGAACTGAGGCCAGAGGCC
CGGCGGTTAGGGGGCCTCCAAATCTGGAGAGCCCTGA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```

Restriction Sites: SgfI-MluI



[View online »](#)

Plasmid Map:



ACCN: NM_024297

Insert Size: 1212 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_024297.2](#)

RefSeq Size: 2041 bp

RefSeq ORF: 1212 bp

Locus ID: 79142

UniProt ID: [Q9BUL5](#)

Cytogenetics: 17p13.1

Protein Families:	Druggable Genome, Transcription Factors
MW:	43.8 kDa
Gene Summary:	Acts as a negative regulator of autophagy, through promoting ubiquitination and degradation of LRSAM1, an E3 ubiquitin ligase that promotes autophagy in response to starvation or infecting bacteria.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (1).