

Product datasheet for **SC317678**

TMEM118 (RNFT2) (NM_032814) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM118 (RNFT2) (NM_032814) Human Untagged Clone
Tag:	Tag Free
Symbol:	TMEM118
Synonyms:	TMEM118
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317678 representing NM_032814. Blue=Insert sequence Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCC**CGATCGCC**
 ATGTGGCTCTTCACAGTGAATCAGGTGTTAAGGAAGATGCAGAGACGCCACAGCAGCAACACGGATAAC
 ATTCCACCTGAAAGAAACCGCAGCCAGGCGCTCAGCTCCGAGGCGAGTGTGGATGAAGGTGGCGTCTTT
 GAGAGTCTGAAGGCAGAGGCAGCCTCCCCACCAGCGCTCTTCTCGGGCTTATCAGGCAGCCTCCCCACC
 AGCTCGTTCCCTCCAGCCTGGTGGCTCCTCGGCTGGCGGCGGGGACGTGTTTCATCCAGATGCC
 GCGTCCAGGGAGGAAGGAGGGGGCGGGGCGAGGGGGGCGCTACCACCACCGCCAGCCCCACCACCAT
 TTCCACCATGGCGGCCACCGGGGGCTCCCTGCTGCAGCACGTGGGTGGGACCACCGGGGGCACTCG
 GAGGAGGGAGGCGACGAGCAGCCTGGGACGCCCGCCCCGCCCTGTCCGAGCTGAAGGCTGTGATCTGC
 TGGCTCCAGAAAGGACTCCCCTTCATCCTGATCCTCCTGGCCAACTGTGCTTTCAGCATAAGCTCGGC
 ATTGCTGTGTGCATCGGGATGGCCAGCACCTTCGCCTATGCCAACTCCACGCTTCGAGAACAGGTCTCA
 CTGAAGGAGAAGAGGTCACTGCTGGTCACTTGTGGATCCTGGCCTTCTGGCGGGGAACACCCTCTAT
 GTGCTTTATACATTACGCTCCAGCAGCTGTACAACAGCCTCATATTCCTGAAGCCCAACCTGGAGATG
 CTGGACTTCTTTGACCTGCTATGGATTGTGGGATCGCAGACTTGTCTGAAGTACATCACCATCGCC
 CTCAAGTGCCATCATCGTGGCCCTGCCAAGATCATCCTGGCTGTCAAGTCCAAGGAAAGTTCTATCTG
 GTCATCGAGGAGCTGAGCCAGCTGTTCCGATCCCTTGTCCCCATCCAGCTGTGGTACAAATACATCATG
 GGTGACGACTCCTCAACAGCTACTTCTGGCGGGGTCTGATCGTTCTCTACAGCCTCTGCAAGTCC
 TTCGACATCTGTGGACGTGTGGCGGAGTTAGGAAAGCCCTGAAGCTTCTCTGTACCTCTCAGAATAT
 GGAGTCCGAGCCACCGGGCAGCAGTGCACAGAAGCTGGTGACATCTGCGCCATCTGTAGGCCGAGTTT
 CGAGAGCCTCTGATTCTCCTGTGCCAGATGCTGTTGAAGGGGCACAAGAAATTGGAGCTGGAGAAGATT
 GATGAAAGTGCAGGTGT**TAA**
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_032814
Insert Size:	1263 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_032814.3</u>
RefSeq Size:	2431 bp
RefSeq ORF:	1263 bp
Locus ID:	84900
UniProt ID:	<u>Q96EX2</u>
Cytogenetics:	12q24.22
Protein Families:	Transmembrane
MW:	46 kDa