

Product datasheet for SC317681

CNPase (CNP) (NM_033133) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	CNPase (CNP) (NM_033133) Human Untagged Clone
Tag:	Tag Free
Symbol:	CNP
Synonyms:	CNP1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317681 representing NM_033133. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGAACAGAGGCTTCTCCGAAAAAGCCACACATTCTGCCCAAGATCTTCTCCGAAGATGTCATCC
TCAGGGGCAAGGACAAGCCTGAGCTGCAGTTTCCCTTCCAGGATGAGGACACAGTGGCCACGCTG
CTAGAGTGCAAGACGCTTTCATCTTGCAGCGCCTGCCAGGAAGCGGCAAGTCCACGCTGGCACGGGTC
ATCGTGGACAAGTACCGTATGGACCAAGATGGTGTCCGCTGACGCTTACAAGATCACCCCGGGCT
CGAGGAGCCTTCCGAGGAGTACAAGCGGCTCGATGAGGACCTGGCTGCCTACTGCCCGCCGGGAC
ATCAGAATCTTGTGCTTGATGACCAACCACGAACGGGAACGGCTGGAGCAGCTTTTGAATGGCC
GACCAGTACCAGTACCAGGTGGTGTGGTGGAGCCCAAGACGGCGTGGCGGCTGGACTGTGCCAGCTC
AAGGAGAAGAACCAGTGGCAGCTGTCCGCTGATGACCTGAAGAAGCTGAAGCCTGGGCTGGAGAAGGAC
TTCTGCCGCTCTACTTCGGCTGGTTTCTGACCAAGAAGACTCTGAGACCCTCCGAAAGCCGGCCAG
GTCTTCTGGAAGAGCTGGGGAACCACAAGGCCTTCAAGAAGGAGCTGCGACAATTCGTCCCTGGGGAT
GAGCCCAGGGAGAAGATGGACTTGGTACCTACTTTGAAAGAGACCCCAAGCGTGTGCATTGCACA
ACCAAGTTTTGTGACTACGGGAAGGCTCCCGGGCAGAGGAGTACGCTCAACAAGATGTGTTAAAGAAA
TCTTACTCCAAGGCCTTACGCTGACCATCTCTGCCCTTTGTGACACCAAGACGACTGGGGCCCGG
GTGGAGTTAAGCGAGCAGCAACTGCAGTTGTGGCCGAGTGTGGACAAGCTGCACCCACTGACCA
CTGCCCGGGGAGCCGCGCCACATCACCTCGGCTGTGACGCTGACGTAGAGGCCGTGACAGCGGC
CTTGACCTCTTAGAGATTCTGCGCAGGAGAAGGGGGCAGCCGAGGCGAGGAGTGGGCGAGCTAAGC
CGGGGCAAGCTCTATTCTTGGGCAATGGGCGCTGGATGCTGACCCTGGCCAAGAACATGGAGGTACGG
GCCATCTTACGGGTACTACGGAAAGGCAAACCTGTGCCACGCAAGGTAGCCGAAGGGGGCGCC
TTGCAGTCTGCACCATCATATGA
ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_033133

Insert Size: 1266 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_033133.4](#)

RefSeq Size: 5222 bp

RefSeq ORF: 1266 bp

Locus ID: 1267

UniProt ID: [P09543](#)

Cytogenetics: 17q21.2

MW: 47.6 kDa

Gene Summary: May participate in RNA metabolism in the myelinating cell, CNP is the third most abundant protein in central nervous system myelin.[UniProtKB/Swiss-Prot Function]