

## Product datasheet for SC317682

### PIP4K2C (NM\_024779) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	PIP4K2C (NM_024779) Human Untagged Clone
Tag:	Tag Free
Symbol:	PIP4K2C
Synonyms:	PIP5K2C
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317682 representing NM_024779. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTGAACCGTCAGAATTTTGTAAACGACTACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGGCGTCCTCCTCGGTCCCACCAGCCACGGTATCGGCGGCGACAGCAGGCCCGGCCAGGTTTCGGC
TTCGCCTCCAAGACCAAGAAGAAGCATTTCGTGCAGCAGAAGGTGAAGGTGTTCCGGGCGGCCGACCCG
CTGGTGGGTGTGTTCTGTGGGGCTAGCCACTCGATCAATGAGCTCAGCCAGGTGCCTCCCCGGTG
ATGCTGCTGCCAGATGACTTTAAGGCCAGCTCCAAGATCAAGGTCAACAATCACCTTTTCCACAGGGAA
AATCTGCCAGTCATTTCAAGTTCAAGGAGTATTGTCCCAGGTCTTCAGGAACCTCCGTGATCGATTT
GGCATTGATGACCAAGATTACTTGGTGTCCCTTACCCGAAACCCCCAGCGAAAGTGAAGGCAGTGAT
GGTCGCTTCCCTTATCTCCTACGATCGGACTCTGGTCATCAAGAAGTATCCAGTGAGGACATTGCTGAC
ATGCATAGCAACCTCTCCAATATCACCAGTACATTGTGAAGTGCCATGGCAACACGCTTCTGCCCCAG
TTCTGGGGATGTACCGAGTCAGTGTGGACAACGAAGACAGCTACATGCTTGTGATGCGCAATATGTTT
AGCCACCGTCTTCTGTGCACAGGAAGTATGACCTCAAGGGTTCCTAGTGTCCCGGGAAGCCAGCGAT
AAGGAAAAGGTTAAAGAATTGCCACCCTTAAGGATATGGACTTTCTCAACAAGAACCAGAAAAGTATAT
ATTGGTGAAGAGGAGAAGAAAATATTTCTGGAGAAGCTGAAGAGAGATGTGGAGTTTCTAGTGCAGCTG
AAGATCATGGACTACAGCCTTCTGCTAGGCATCCACGACATCATTCGGGGCTCTGAACCAGAGGAGGAA
GCGCCCGTGGGAGGATGAGTCAGAGGTGGATGGGGACTGCAGCCTGACTGGACCTCCTGCTCTGGTG
GGCTCCTATGGCACCTCCCAGAGGGTATCGGAGGCTACATCCATTCCCATCGGCCCTGGGCCAGGA
GAGTTTGAAGTCTTCAATTGATGTCTATGCCATCCGGAGTGTGAAGGAGCCCCCAGAAGGAGGTCTAC
TTCATGGGCTCATTGATATCCTTACACAGTATGATGCCAAGAAGAAAGCAGCTCATGCAGCCAAAAT
GTCAAGCATGGGGCTGGGCAGAGATCTCTACTGTCCATCCGGAGCAGTATGCTAAGCGATTCTGGAT
TTTATTACCAACATCTTTGCC TAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

**Plasmid Map:**



ACCN: NM\_024779

Insert Size: 1266 bp

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

**OTI Annotation:** This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_024779.4](#)

RefSeq Size: 3229 bp

RefSeq ORF: 1266 bp

Locus ID: 79837

UniProt ID:	<a href="#">Q8TBX8</a>
Cytogenetics:	12q13.3
Domains:	PIP5K
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism, Phosphatidylinositol signaling system, Regulation of actin cytoskeleton
MW:	47.3 kDa
Gene Summary:	May play an important role in the production of Phosphatidylinositol bisphosphate (PIP2), in the endoplasmic reticulum.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (1) represents the longest transcript. Both variants 1 and 2 encode the same isoform (a).