

Product datasheet for **SC317683**

UBE2Q1 (NM_017582) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	UBE2Q1 (NM_017582) Human Untagged Clone
Tag:	Tag Free
Symbol:	UBE2Q1
Synonyms:	GTAP; NICE-5; NICE5; PRO3094; UBE2Q
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317683 representing NM_017582. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCCGCATCGCC
ATGCAGCAGCCGAGCCGAGGGGCGAGCAGCAGCCGGGGCCGGGGCAGCAGCTGGGGGGCCAGGGGGCG
GCGCCGGGGGGCCGGGGCGGCCAGGGGGGGCCGGGGCCGGGGCCCTGCCTGAGGCGAGAGCTGAAG
CTGCTCGAGTCCATCTTCCACCGCGGCCACGAGCGCTTCCGCATTGCCAGCGCTGCCTGGACGAGCTG
AGCTGCGAGTTCTGCTGGCTGGGGCCGAGGGGCGGGGCGGGGCGCCCGCCCGACCGCATCTCCCC
CCACGGGGGTCGGTGCCTGGGGATCCTGTCCGCATCCACTGCAACATCACGGAGTCATACCCTGCTGTG
CCCCCATCTGGTGGTGGAGTCTGATGACCCTAACTGGCTGCTGTCTTGGAGAGGCTGGTGGACATA
AAGAAAGGAATACTCTGCTATTGAGCATCTGAAGAGGATCATCTCCGACCTGTGTAAACTCTATAAC
CTCCCTCAGCATCCAGATGTGGAGATGCTGGATCAACCCTTGCCAGCAGAGCAGTGCACACAGGAAGAC
GTGTCTTCAGAAGATGAAGATGAGGAGATGCCTGAGGACACAGAAGACTTAGATCACTATGAAATGAAA
GAGGAAGAGCCAGCTGAGGGCAAGAAATCTGAAGATGATGGCATTGGAAGAAAGAAACTTGGCCATCCTA
GAGAAAATTAAGAAGAACAGAGGCAAGATTACTTAAATGGTGCAGTGTCTGGCTCGGTGCAGGCCACT
GACCGGCTGATGAAGGAGCTCAGGGATATATACCGATCACAGAGTTTCAAAGCGGAAACTATGCAGTC
GAACTCGTGAATGACAGTCTGTATGATTGGAATGTCAAACCTCAAAGTTGACCAGGACAGCGCTTTG
CACACGATCTCCAGATCCTCAAAGAGAAAGAAGGAGCCGACTTACTTAACTTTCTTTAAA
GATAACTTTCCCTTTGACCCACCATTTGTCAGGGTTGTGTCTCCAGTCTCTCTGGAGGGTATGTTCTG
GGCGGAGGGGCCATCTGCATGAACTTCTCACAAACAGGGCTGGAGCAGTGCCTACTCCATAGAGTCA
GTGATCATGCAGATCAGTGCCACTGGTGAAGGGGAAAGCACGAGTGCAGTTTGGAGCCAAACAACTCT
CAATACAGTCTGACAAGAGCACAGCAGTCTACAAGTCTTGGTGCAGATCCACGAAAAAACCGCTGG
TACACACCCCAAAAGAAGACGGCTAA
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online >](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_017582

Insert Size: 1269 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_017582.6](#)

RefSeq Size: 3223 bp

RefSeq ORF: 1269 bp

Locus ID: 55585

UniProt ID: [Q7Z7E8](#)

Cytogenetics: 1q21.3

Domains: UBCc

Protein Pathways: Ubiquitin mediated proteolysis

MW: 46.1 kDa

Gene Summary: The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes (E1s), ubiquitin-conjugating enzymes (E2s), and ubiquitin-protein ligases (E3s). This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. The encoded protein is 98% identical to the mouse counterpart. [provided by RefSeq, Jul 2008]