

## Product datasheet for **SC317690**

### DIPK1A (NM\_001006605) Human Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	DIPK1A (NM_001006605) Human Untagged Clone
Tag:	Tag Free
Symbol:	DIPK1A
Synonyms:	FAM69A
Mammalian Cell Selection:	None
Vector:	<u>pCMV6-XL4</u>
E. coli Selection:	Ampicillin (100 ug/mL)
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001006605, the custom clone sequence may differ by one or more nucleotides

```

ATGGCGAGGAGTCTCTGTCGGGGCCCTGGCTAAGGAAACCTATTACCTCCAGGCTCGC
TTCTCATATGTGCGGATGAAATATCTTTTCTTTCTGGTTAGTGTTTTGTTGGAAGC
TGGATTATATATGTGCACTATTCTACCTATACAGAATTATGCAGAGGAAAGGACTGTAAG
AAAATAATATGTGACAAGTACAAGACTGGAGTTATTGATGGGCCTGCATGTAACAGCCTT
TGTGTTACAGAACTCTTTACTTTGGAAAATGTTTATCCACCAAGCCCAACAATCAGATG
TATTTAGGGATTTGGGATAATCTACCAGGTGTTGTGAAATGTCAAATGGAACAAGCGCTT
CATCTTGATTTTGGAACTGAATTGGAACCAAGAAAAGAAATAGTGCTATTTGATAAGCCA
ACTAGAGGAACTACTGTACAAAAATTTAAAGAAATGGTCTATAGTCTCTTTAAGGCAAAA
TTGGGTGACCAAGGAAACCTCTCTGAAGTGGTTAATCTCATCTTGACGGTGGCTGATGGA
GACAAAGATGGCCAGGTTTCCTTGGGAGAAGCAAAGTCGGCATGGGCACTTCTTCAACTG
AATGAATTTCTTCTCATGGTGATACTTCAAGATAAAGAACATACCCCAAATTAATGGGA
TTCTGTGGTGACCTCTATGTGATGGAAAGTGTGAATATACCTCTCTTTATGGAATAAGC
CTTCTTGGGTCATTGAACTTTTTATTCCATCTGGGTTCCAGAAGAAGCATGGATCAGCTG
TTCACACCATCATGGCCAAGAAAGGCCAAAATAGCCATAGGACTTCTAGAATTTGTGGAA
GGATATAATGATAAGTATGATTTGAAAAATGGTGGATATGAGAAAAATTGTCCAGAGACA
AACCTGAAAGAACTTATTAAGGATCGTCACTGTGAGTCTGATTTGGACTGTGCTATGGC
ACAGATTGTAGAACTAGCTGTGATCAGAGTACAATGAAGTGTACTTCAGAAGTGATACAA
CCAAACTTGCCAAAAGCTTGTCAAGTACTCAAAGACTACCTACTGCGTGGTGCTCCAAGT
GAAATTCGTGAAGAATTAGAAAAGCAGCTTTATTCTTGTATTGCTCTCAAAGTCACAGCA
AATCAAATGGAATGGAACATTCTTTGATACTAAATAACCTAAAAACATTATTGTGGAAG
AAAATTTCTACACTAATGACTCT

```

Restriction Sites: Please inquire



[View online »](#)

<b>ACCN:</b>	NM_001006605
<b>Insert Size:</b>	2603 bp
<b>OTI Disclaimer:</b>	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at <a href="mailto:custsupport@origene.com">custsupport@origene.com</a> or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a></p>
<b>OTI Annotation:</b>	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u><a href="#">NM_001006605.3</a></u> , <u><a href="#">NP_001006606.2</a></u>
<b>RefSeq Size:</b>	2603 bp
<b>RefSeq ORF:</b>	1287 bp
<b>Locus ID:</b>	388650
<b>UniProt ID:</b>	<u><a href="#">Q5T7M9</a></u>
<b>Cytogenetics:</b>	1p22.1
<b>Protein Families:</b>	Transmembrane

**Gene Summary:**

This gene encodes a member of the FAM69 family of cysteine-rich type II transmembrane proteins. These proteins localize to the endoplasmic reticulum but their specific functions are unknown. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Nov 2011]

Transcript Variant: This variant (1) encodes the longest isoform (1). Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.