

Product datasheet for SC317701

ZNF829 (NM_001037232) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	ZNF829 (NM_001037232) Human Untagged Clone
Tag:	Tag Free
Symbol:	ZNF829
Vector:	<u>pCMV6 series</u>
Fully Sequenced ORF:	>NCBI ORF sequence for NM_001037232, the custom clone sequence may differ by one or more nucleotides

```

ATGCCCCATTCTCCTCTGATCTCCATTCTCATGTGTGGTGTACCCAGAAGAGGAGGAA
AGAATGCATGATGAATTCTACAAGCAGTATCCAAGGGGCCGGTGATGTTCAAGGATGTT
TCCATAGACTTCTCTCAAGAGGAATGGGAATGCCTGGACGCTGATCAGATGAATTATAC
AAAGAAGTGATGTTGGAGAATTTCAACACCTGGTTTCAGTGGGACTTTCCAATTCTAAG
CCAGCTGTGATCTCCTTATTGGAACAAGGAAAAGAGCCCTGGATGGTTGATAGAGAGCTG
ACTAGAGGCCTGTGTTCAAGTCAATGTGTGAGACCAAAATATTATCTCTAAAG
AAGAGACATTTTCAGTCAAGTAATAATTACCCGTGAAGACATGTCTACTTTTATTCAGCCC
ACATTTCTTATTCCACCTCAAAAACTATGAGTGAAGAGAAACCATGGGAATGTAAGATA
TGTGGAAGACCTTTAATCAAACTCACAATTTATCCAACATCAGAGAATTCATTTTGGT
GAAAAACACTATGAATCTAAGGAGTATGGGAAGTCCTTTAGTCGTGGCTCACTCGTTACT
CGACATCAGAGGATTCACACTGGTAAAAAACCTATGAATGTAAGGAATGTGGCAAGGCT
TTTAGTTGTAGTTCATATTTTCTCAACATCAGAGGATTCACACTGGTGAGAAACCTAT
GAATGTAAGGAATGTGGAAAAGCCTTTAAGTATTGCTCAAACCTAATGATCATCAGAGA
ATTCACACTGGTGAGAAACCTATGAATGTAAGTATGTGGAAAAGCCTTTACTAAAAGT
TCACAATTTTTCTACATCTGAGAATTCATACTGGTGAGAAACCTATGAATGTAAGAA
TGTGGGAAAGCCTTTACTCAACACTCAAGGCTTATTCAGCATCAGAGAATGCATACTGGT
GAGAAACCTTATGAATGTAAGCAGTGTGGGAAGGCCTTTAATAGTGCCTCAACACTTACT
AACCATCAGAGAATTCATGCTGGTGAGAAGCTCTATGAATGTGAAGAATGTAGAAAGGCC
TTTATTCAGAGCTCAGAACTTATTCAACATCAGAGAATCCATACAGATGAAAAACCATAT
GAATGTAATGAATGTGGGAAGGCCTTTAATAAGGCTCAAATCTTACTCGACATCAGAGA
ATTCACACTGGTGAGAAACCTATGACTGTAAGGAATGTGGAAAGGCTTTTGGTAGTCGC
TCTGACCTCATTCGCCATGAGGAATTCATACTGGT
  
```

Restriction Sites:	Please inquire
ACCN:	NM_001037232
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037232.2</u> , <u>NP_001032309.2</u>
RefSeq Size:	1802 bp
RefSeq ORF:	1299 bp
Locus ID:	374899
UniProt ID:	<u>Q3KNS6</u>
Cytogenetics:	19q13.12
Gene Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function] Transcript Variant: This variant (2) differs in the 5' UTR, lacks a portion of the 5' coding region, and uses a downstream start codon, compared to variant 1. The encoded isoform (2) is shorter at the N-terminus, compared to isoform 1. Sequence Note: This RefSeq record was created from transcript and genomic sequence data to make the sequence consistent with the reference genome assembly. The genomic coordinates used for the transcript record were based on transcript alignments.