

Product datasheet for SC317711

GCET1 (SERPINA9) (NM_175739) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	GCET1 (SERPINA9) (NM_175739) Human Untagged Clone
Tag:	Tag Free
Symbol:	SERPINA9
Synonyms:	GCET1; SERPINA11; SERPINA11b
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317711 representing NM_175739. Blue=Insert sequence Red=Cloning site Green=Tag(s)


```
GCTCGTTTGTAGTAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGCCGGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCAAGGACAGGGCAGGAGAAGAGGAACCTGCAAAGACATATTTTGTCCAAAATGGCATCTTACCTT
TATGGAGTACTCTTTGCTGTTGGCCTCTGTGCTCCAATCTACTGTGTGTCCCCGGCCAATGCCCCAGT
GCATACCCCGCCCTTCTCCACAAGAGCACCCCTGCCTCACAGGTGTATTCCCTCAACACCGACTTT
GCCTTCCGCCTATACCGCAGGCTGGTTTTGGAGACCCCGAGTCAGAACATCTTCTCTCCCTGTGAGT
GTCTCCACTTCCCTGGCCATGCTCTCCCTGGGGCCCACTCAGTCACCAAGACCCAGATTCTCCAGGGC
CTGGGCTTCAACCTCACACACACCAGAGTCTGCCATCCACCAGGGCTCCAGCACCTGGTTCACTCA
CTGACTGTTCCCAGCAAAGACCTGACCTTGAAGATGGGAAGTGCCTCTTCGTCAGAAGGAGCTGCAG
CTGCAGGCAAATTTCTTGGGCAATGTCAAGAGGCTGTATGAAGCAGAAGTCTTTTCTACAGATTTCTCC
AACCCCTCCATTGCCAGGCGAGGATCAACAGCCATGTGAAAAAGAACCCAAAGGGAAGGTTGTAGAC
ATAATCCAAGGCCTTGACCTTCTGACGGCCATGGTTCTGGTGAACCACATTTTCTTTAAAGCCAAGTGG
GAGAAGCCCTTACCCTGAATATACAAGAAAGAACTCCCATCCTGGTGGGCGAGCAGGTCACTGTG
CATGTCCCATGATGCACCAGAAAGAGCAGTTCGCTTTGGGGTGGATACAGAGCTGAAGTCTTTGTG
CTGCAGATGGATTACAAGGGAGATGCCGTGGCCTCTTTGTCTCCCTAGCAAGGCAAGATGAGGCAA
CTGGAACAGGCCTTGTGAGCCAGAACACTGAGAAAGTGGAGCCACTCACTCCAGAAAAGGTGATAGAG
GTGTTTCATCCCCAGATTTCCATTTCTGCCTCCTACAATCTGAAACCCTCCTCCCGAAGATGGGCATC
CAAAATGTCTTTGACAAAATGCTGATTTTTCTGGAATTGCAAAGAGAGACTCCCTGCAGGTTTCTAAA
GCAACCCACAAGGCTGTGCTGGATGTCAGTGAAGAGGGCACTGAGGCCACAGCAGCTACCACCACCAAG
TTCATAGTCCGATCGAAGGATGGCCCTTACTTCACTGTCTCCTTCAATAGGACCTTCTGATGATG
ATTACAATAAAGCCACAGACGGTATTCTTTCTAGGAAAAGTGGAAAATCCCCTAAATCCTAG
ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)

Restriction Sites: SgfI-MluI

Plasmid Map:



ACCN: NM_175739

Insert Size: 1308 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

OTI Annotation: This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175739.3](#)

RefSeq Size: 1851 bp

RefSeq ORF: 1308 bp

Locus ID: 327657

UniProt ID: [Q86WD7](#)

Cytogenetics: 14q32.13

Protein Families: Druggable Genome, Secreted Protein, Transmembrane

MW: 48.6 kDa

Gene Summary: Protease inhibitor that inhibits trypsin and trypsin-like serine proteases (in vitro). Inhibits plasmin and thrombin with lower efficiency (in vitro).[UniProtKB/Swiss-Prot Function]
Transcript Variant: This variant (1) represents the longest transcript and encodes the longest isoform (A).