

Product datasheet for SC317712

PTD016 (RNFT1) (NM_016125) Human Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PTD016 (RNFT1) (NM_016125) Human Untagged Clone
Tag:	Tag Free
Symbol:	PTD016
Synonyms:	PTD016
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>SC317712 representing NM_016125. Blue=Insert sequence Red=Cloning site Green=Tag(s)

```

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
ATGCCGCTGTTCTTGCTGTCGCTCCCGACACCTCCGTCGCTTCTGGTCATGAGAGGAGACAGAGGCCT
GAAGCAAAGACATCTGGGTGAGAGAAAAGTATTTAAGGGCCATGCAAGCCAATCGTAGCCAACCTGCAC
AGTCCTCCAGGAAGTGAAGCAGTGAGGATGCCTCAACCCCTCAGTGTGTCCACACAAGATTGACAGGA
GAGGGTCTTGCCCTCATTCTGGAGATGTTTATCCAGATAAACTCCATACCTAAAGAATGTGCAGAA
AATGCAAGCTCCAGAAATATAAGGTCAGGTGTCCATAGCTGTGCCATGGATGTGTACACAGTCGCTTA
CGGGGTCACTCCACAGTGAAGCAAGGCTGACTGATGATACTGCCGAGAATCTGGAGATCATGGTAGT
AGCTCCTTCTCAGAATTCGCTATCTCTCAAGTGGCTGCAAAAAGTCTTCATATATTTTGATTCTG
AGCGTCAAACCTGTTATGCAGCATATAACAGGAATTTCTCTTGAATTTGGGCTGCTAACAACTTTTATG
TATGCAACAAAAGCATTGTAAATCAGGTTTTCTAAGAGAAAGGTCCTCAAAGATTCAAGTGTGCTTGG
TTACTGGTATTCTTAGCAGGATCTTCTGTTCTTTTATATTACACCTTTCATTCTCAGTCACTTTATTAC
AGCTTAATTTTTTAAATCCTACTTTGGACCATTTGAGCTTCTGGGAAGTATTTTGGATTGTTGGAATT
ACAGAGTTCATTCTGAAATCTTTTTTATGGGCTTAAATGCCTTATTTTATTGGTGCCTTCTTTTATC
ATGCCCTTTTAAATCTAAGGTTACTGGTATATGCTTTTAGAAGAATTGTGTCAATACTACCGAATTTT
GTTCCCATACAGTTTGGTTTCGCTACCTTATAAGCTATGGGGAGTTTGGTAACGTAACCTAGATGGAGT
CTTGGGATACTGCTGGCTTTACTCTACCTCATATTTAAACTTTTGGAAATTTTGGGCATCTGAGAACT
TTCAGACAGGTTTACGAATATTTTACACAACCAAGTTATGGAGTGGCTGCCAGCAAGAGACAGTGT
TCAGATGTGGATGATATTTGTTCAATATGTCAAGCTGAATTTGAGAAGCAATCTTCTCATTTGTCAG
CATATATTTTGTGAAGAGTGCATGACCTTATGGTTTAAACAGAGAGAAAACATGTCCACTCTGCAGAACT
GTGATTTGAGACCATATAAACAATGGAAGGATGGAGCCACTTCATCACACCTTCAAATATATTAA
ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
  
```


[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_016125
Insert Size:	1308 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
OTI Annotation:	This TrueClone is provided through our Custom Cloning Process that includes sub-cloning into OriGene's pCMV6 vector and full sequencing to provide a non-variant match to the expected reference without frameshifts, and is delivered as lyophilized plasmid DNA.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_016125.3</u>
RefSeq Size:	2031 bp
RefSeq ORF:	1308 bp
Locus ID:	51136
UniProt ID:	<u>Q5M7Z0</u>
Cytogenetics:	17q23.1
Protein Families:	Druggable Genome, Transmembrane
MW:	49.7 kDa
Gene Summary:	E3 ubiquitin-protein ligase that acts in the endoplasmic reticulum (ER)-associated degradation (ERAD) pathway, which targets misfolded proteins that accumulate in the endoplasmic reticulum (ER) for ubiquitination and subsequent proteasome-mediated degradation. Protects cells from ER stress-induced apoptosis.[UniProtKB/Swiss-Prot Function]